

Product Description SALSA® MLPA® Probemix P197-A4 KCNQ3

To be used with the MLPA General Protocol.

Version A4. Compared to version A3, two reference probes have been replaced and a small change in length of two target probes. For complete product history see page 7.

Catalogue numbers:

- **P197-025R:** SALSA MLPA Probemix P197 KCNQ3, 25 reactions.
- **P197-050R:** SALSA MLPA Probemix P197 KCNQ3, 50 reactions.
- **P197-100R:** SALSA MLPA Probemix P197 KCNQ3, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mlpa.com).

Certificate of Analysis: Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mlpa.com.

Precautions and warnings: For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mlpa.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information: The SALSA MLPA Probemix P197 KCNQ3 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *KCNQ3*, *CHRNA4*, *EPM2A*, *NHLRC1* (also known as *EPM2B*), and *CHRN2* genes, which are associated with Epilepsy.

Defects in the *KCNQ3* gene can cause benign familial neonatal convulsions type 2 (BFNC2) or benign neonatal type 2 (EBN2), which is also known as epilepsy. The *KCNQ3* gene comprises 15 exons and spans ~360 kb of genomic DNA on chromosome 8q24.22. The *CHRNA4* gene comprises 6 exons spanning ~18 kb of genomic DNA on chromosome 20q13.33. Mutations in this gene appear to account for a small proportion of the cases of nocturnal frontal lobe epilepsy. Mutations in the *EPM2A* gene have been associated with myoclonic epilepsy of Lafora. The *EPM2A* gene encodes the laforin protein, comprises 4 exons, and spans ~111 kb of genomic DNA on chromosome 6q24.3. The *NHLRC1* (*EPM2B*) gene encoding the malin protein, comprises 1 exon and spans ~1.2 kb of genomic DNA on chromosome 6p22.3. Mutations in *NHLRC1* cause progressive myoclonus epilepsy and are also associated with autosomal dominant nocturnal frontal lobe epilepsy. The *CHRN2* gene comprises 6 exons spanning ~12.2 kb of genomic DNA on chromosome 1q21.3. Finally, the *KCNQ1* gene comprises 17 exons and spans ~404 kb of genomic DNA on chromosome 11p15.5-4.

This SALSA MLPA Probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>

For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>

Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering: The *CHRN2*, *NHLRC1*, *EPM2A*, *KCNQ3*, *KCNQ1*, and *CHRNA4* exon numberings used in this P197-A4 KCNQ3 product description are the exon numberings from the RefSeq transcripts NM_000748.3, NM_198586.3, NM_005670.4, NM_004519.4, NM_000218.3, and NM_000744.6, which are identical to the NG_008027.1, NG_016750.1, NG_012832.2, NG_008854.2, LRG_287, and NG_011931.1 sequences respectively. The exon numbering and NM_ sequence used have been retrieved on 02/2020. As changes to the NCBI database can occur after release of this product description, exon numbering may not be up-to-date.

Probemix content: The SALSA MLPA Probemix P197-A4 KCNQ3 contains 39 MLPA probes with amplification products between 137 and 463 nucleotides (nt). This includes 15 probes for the *KCNQ3* gene, six probes for the *CHRNA4* gene, four probes for the *EPM2A* gene, two probes for *NHLRC1* gene, two probes for the *CHRNA2* gene, and finally two probes for the *KCNQ1* gene. In addition, eight reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mlpa.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mlpa.com.

| Length (nt) | Name |
|-------------|--|
| 64-70-76-82 | Q-fragments (only visible with <100 ng sample DNA) |
| 88-96 | D-fragments (low signal of 88 nt and 96 nt fragment indicates incomplete denaturation) |
| 92 | Benchmark fragment |
| 100 | X-fragment (X chromosome specific) |
| 105 | Y-fragment (Y chromosome specific) |

MLPA technique: The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mlpa.com).

MLPA technique validation: Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens: Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples: A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from unrelated individuals who are from families without a history of epilepsy. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol.

Positive control DNA samples: MRC-Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de/home.html>) have a diverse collection of biological resources which may be used as a positive control DNA sample in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis: Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mlpa.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results: The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the dosage quotient (DQ) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the DQ of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

| Copy number status | Dosage quotient |
|--|--------------------|
| Normal | $0.80 < DQ < 1.20$ |
| Homozygous deletion | $DQ = 0$ |
| Heterozygous deletion | $0.40 < DQ < 0.65$ |
| Heterozygous duplication | $1.30 < DQ < 1.65$ |
| Heterozygous triplication/Homozygous duplication | $1.75 < DQ < 2.15$ |
| Ambiguous copy number | All other values |

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Incomplete DNA denaturation (e.g. due to salt contamination) can lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- When running MLPA products, the capillary electrophoresis protocol may need optimization. False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: lower injection voltage / injection time settings, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure:

- In most populations, the major cause of genetic defects in the *KCNQ3*, *CHRNA4*, *EPM2A*, *NHLRC1*, *KCNQ1* and *CHRNA2* genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P197 KCNQ3.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region do exist but remain undetected.
- Sequence changes (e.g. SNPs, point mutations, small indels) in the target sequence detected by a probe can cause false positive results. Mutations/SNPs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results: Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

LOVD mutation database: <https://databases.lovd.nl/shared/genes>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNPs and unusual results (e.g., a duplication of *KCNQ3* exons 4 and 6 but not exon 5) to MRC-Holland: info@mlpa.com.

Table 1. SALSA MLPA Probemix P197-A4 KCNQ3

| Length (nt) | SALSA MLPA probe | Chromosomal position (hg18) ^a | | | |
|-------------|--|--|-----------------|-----------------|----------------|
| | | Reference | KCNQ3 | EPM2A | CHRNA4 |
| 64-105 | Control fragments – see table in probemix content section for more information | | | | |
| 137 | Reference probe 03797-L04594 | 21q22 | | | |
| 142 | KCNQ3 probe 06595-L06153 | | Exon 3 | | |
| 154 | Reference probe 20337-L27719 | 1p36 | | | |
| 160 Δ | EPM2A probe 06617-L29224 | | | Exon 1 | |
| 166 | KCNQ3 probe 06596-L06154 | | Exon 4 | | |
| 172 | KCNQ3 probe 06603-L06161 | | Exon 11 | | |
| 178 | CHRN2 probe 06616-L06174 | | | | Exon 6 |
| 184 | KCNQ3 probe 06601-L06159 | | Exon 9 | | |
| 190 * | Reference probe 22509-L31658 | 14q32 | | | |
| 198 | KCNQ3 probe 06604-L29191 | | Exon 12 | | |
| 205 | NHLRC1 probe 06621-L07195 | | | | Exon 1 |
| 214 | CHRNA4 probe 06609-L06167 | | | Exon 2 | |
| 220 | CHRNA4 probe 06613-L06171 | | | Exon 6 | |
| 228 † | EPM2A probe 06619-L32046 | | | Exon 3 | |
| 234 * | Reference probe 11156-L16377 | 5q31 | | | |
| 240 † | CHRNA4 probe 06610-L10371 | | | Exon 3 | |
| 247 | EPM2A probe 06618-L06176 | | | Exon 2 | |
| 256 | EPM2A probe 06620-L07196 | | | Exon 4 | |
| 265 | Reference probe 03241-L02678 | 13q14 | | | |
| 274 ‹ | KCNQ3 probe 06593-L06151 | | Upstream | | |
| 283 | KCNQ3 probe 06600-L29155 | | Exon 8 | | |
| 293 ‹ | CHRN2 probe 06615-L06173 | | | | Exon 2 |
| 301 | KCNQ1 probe 03551-L02917 | | | | Exon 14 |
| 310 | KCNQ3 probe 06597-L06155 | | Exon 5 | | |
| 320 | KCNQ3 probe 06602-L06160 | | Exon 10 | | |
| 328 | NHLRC1 probe 06622-L06180 | | | | Exon 1 |
| 337 | Reference probe 04097-L02899 | 7q36 | | | |
| 346 | CHRNA4 probe 06611-L07197 | | | Exon 4 | |
| 364 | CHRNA4 probe 06612-L06170 | | | Exon 5 | |
| 378 | Reference probe 05921-L05366 | 17q11 | | | |
| 388 | KCNQ3 probe 06594-L29158 | | Exon 2 | | |
| 394 | KCNQ3 probe 06606-L29192 | | Exon 14 | | |
| 400 | CHRNA4 probe 06608-L06166 | | | Upstream | |
| 409 | KCNQ1 probe 03555-L02921 | | | | Exon 17 |
| 418 | KCNQ3 probe 06598-L06156 | | Exon 6 | | |
| 427 | KCNQ3 probe 06607-L07198 | | Exon 15 | | |
| 436 | KCNQ3 probe 07315-L06163 | | Exon 13 | | |
| 454 ‹ | KCNQ3 probe 08195-L08089 | | Exon 1 | | |
| 463 | Reference probe 11713-L12484 | 10q22 | | | |

a) See above section on exon numbering for more information.

* New in version A4.

† Changed in version A4. Minor alteration, no change in sequence detected.

‹ Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Table 2. P197-A4 probes arranged according to chromosomal location

Table 2a. *CHRN2*

| Length (nt) | SALSA MLPA probe | <i>CHRN2</i> exon ^a | Ligation site NM_000748.3 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|--------------------------------|---------------------------|---|------------------------|
| | | <i>start codon</i> | 268-270 (Exon 1) | | |
| 293 « | 06615-L06173 | Exon 2 | 440-441 | GCTGGTGACAGT-ACAGCTTATGGT | 6.4 kb |
| 178 | 06616-L06174 | Exon 6 | 1844-1845 | TTGGGTGGAGGA-TGGACGAGTGAG | |
| | | <i>stop codon</i> | 1774-1776 (Exon 6) | | |

Table 2b. *NHLRC1*

| Length (nt) | SALSA MLPA probe | <i>NHLRC1</i> exon ^a | Ligation site NM_198586.3 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|---------------------------------|---------------------------|---|------------------------|
| | | <i>start codon</i> | 72-74 (Exon 1) | | |
| 205 | 06621-L07195 | Exon 1 | 712-713 | CATTGGAGGCCA-ATTCTCCTTACC | 0.4 kb |
| 328 | 06622-L06180 | Exon 1 | 1065-1066 | ATCACCAGGGAA-ATGTGATTGTTG | |
| | | <i>stop codon</i> | 1257-1259 (Exon 1) | | |

Table 2c. *EPM2A*

| Length (nt) | SALSA MLPA probe | <i>EPM2A</i> exon ^a | Ligation site NM_005670.4 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|--------------------------------|---------------------------|---|------------------------|
| | | <i>start codon</i> | 23-25 (Exon 1) | | |
| 160 Δ | 06617-L29224 | Exon 1 | 121 nt before exon 1 | GATGCATCCCAA-AGAAGGCCGAGA | 49.4 kb |
| 247 | 06618-L06176 | Exon 2 | 388-389 | GATGGTGTGTAT-TGTCTCCCAATA | 50.9 kb |
| 228 | 06619-L32046 | Exon 3 | 644-645 | GCTGTAACCGCT-ACCCAGAGCCCA | 7.8 kb |
| 256 | 06620-L07196 | Exon 4 | 942-943 | CTACATTGACGA-AGAGGCCTTGGC | |
| | | <i>stop codon</i> | 1016-1018 (Exon 4) | | |

Table 2d. *KCNQ3*

| Length (nt) | SALSA MLPA probe | <i>KCNQ3</i> exon ^a | Ligation site NM_004519.4 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|--------------------------------|---------------------------|---|------------------------|
| | | <i>start codon</i> | 564-566 (Exon 1) | | |
| 274 « | 06593-L06151 | Upstream (Exon 1) | 961 nt before exon 1 | TTGGAGGGCTCT-CTGGACATTTAC | 1.8 kb |
| 454 « | 08195-L08089 | Exon 1 | 871-872 | AAACAACGCCAA-GTACCGCGCAT | 294.1 kb |
| 388 | 06594-L29158 | Exon 2 | 973-974 | CCTGGGGTGCTT-GATTCTGGCTGT | 1.8 kb |
| 142 | 06595-L06153 | Exon 3 | 1104-1105 | CTGCTGGATGTT-GCTGCCGATACA | 4.1 kb |
| 166 | 06596-L06154 | Exon 4 | 1280-1281 | CGCATGCTGCGG-ATGGACCGGAGA | 4.7 kb |
| 310 | 06597-L06155 | Exon 5 | 1438-1439 | GGTGGATGCACA-AGGAGAGGAGAT | 1.2 kb |
| 418 | 06598-L06156 | Exon 6 | 1497-1498 | TCTCTTTCAGA-TCACACTGGCCA | 3.9 kb |
| 283 | 06600-L29155 | Exon 8 | 1719-1720 | CCTGGAGGTATT-ATGCTACCAACC | 7.0 kb |
| 184 | 06601-L06159 | Exon 9 | 15 nt after exon 9 | AGTTTCTGATTA-TGAATCCCTTC | 22.3 kb |
| 320 | 06602-L06160 | Exon 10 | 2026-2027 | GCAGAGTTCTGA-AGGTAATGCCTT | 1.0 kb |
| 172 | 06603-L06161 | Exon 11 | 2074-2075 | GGGCTATGGGAA-TGACTTCCCCAT | 2.2 kb |
| 198 | 06604-L29191 | Exon 12 | 2191-2192 | GCCTTACGATGT-GAAGGATGTGAT | 3.6 kb |
| 436 | 07315-L06163 | Exon 13 | 2334-2335 | AAGGGTCAGCAT-TCACCTTCCCAT | 2.1 kb |
| 394 | 06606-L29192 | Exon 14 | 2393-2394 | AGACCATCCACA-TCAGAAATCGAA | 2.4 kb |
| 427 | 06607-L07198 | Exon 15 | 2622-2623 | TCATCTGCAACT-ATTCTGAGACAG | |
| | | <i>stop codon</i> | 3180-3182 (Exon 15) | | |

Table 2e. *KCNQ1*

| Length (nt) | SALSA MLPA probe | <i>KCNQ1</i> exon ^a | Ligation site NM_000218.3 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|--------------------------------|---------------------------|---|------------------------|
| | | <i>start codon</i> | 92-94 (Exon 1) | | |
| 301 | 03551-L02917 | Exon 14 (13) | 1749-1750 | CCTCAACCTCAT-GGTGCGCATCAA | 72.8 kb |
| 409 | 03555-L02921 | Exon 17 (16) | 2908-2909 | CCAAACACACAG-AAGGGGACTGCC | |
| | | <i>stop codon</i> | 2120-2122 (Exon 17) | | |

Table 2f. *CHRNA4*

| Length (nt) | SALSA MLPA probe | <i>CHRNA4</i> exon ^a | Ligation site NM_00744.6 | Partial sequence ^b (24 nt adjacent to ligation site) | Distance to next probe |
|-------------|------------------|---------------------------------|--------------------------|---|------------------------|
| | | <i>start codon</i> | 232-234 (Exon 1) | | |
| 400 | 06608-L06166 | Upstream (Exon 1) | 666 nt before exon 1 | CTGCACACGAGA-TTCAGCCGCACA | 2.4 kb |
| 214 | 06609-L06167 | Exon 2 | 364-365 | TGAAGAACTCT-TCTCCGTTACA | 3.3 kb |
| 240 | 06610-L10371 | Exon 3 | 486-487 | ATGATGACCACG-AACGTATGGGTG | 0.4 kb |
| 346 | 06611-L07197 | Exon 4 | 562-563 | ATGTCACCTCCA-TCCGCATCCCCT | 5.2 kb |
| 364 | 06612-L06170 | Exon 5 | 794-795 | CGACAAGGCCAA-GATCGACCTGGT | 4.3 kb |
| 220 | 06613-L06171 | Exon 6 | 2314-2315 | TGTGGAGCTGCT-TCCAGTTGGACT | |
| | | <i>stop codon</i> | 2113-2115 (Exon 6) | | |

a) See above section on exon numbering for more information.

b) Only partial probe sequences are shown. Complete probe sequences are available at www.mlpa.com. Please notify us of any mistakes: info@mlpa.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Related SALSA MLPA probemixes

- P114 Long-QT Contains probes for the *KCNQ1* gene.
- P137 SCN1A Contains probes for the *SCN1A* gene, involved in epilepsy.
- P166 KCNQ2 Contains probes for the *KCNQ2* gene, involved in benign familial neonatal convulsion (BFNC).

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P197 KCNQ3

- Grinton et al. (2015). Familial neonatal seizures in 36 families: Clinical and genetic features correlate with outcome. *Epilepsia* 56:1071-80.
- Soldovieri et al. (2014). Novel KCNQ2 and KCNQ3 Mutations in a Large Cohort of Families with Benign Neonatal Epilepsy: First Evidence for an Altered Channel Regulation by Syntaxin-1A. *Hum Mut.* 35:356-67.

| P197 Product history | |
|----------------------|--|
| Version | Modification |
| A4 | Two reference probes have been replaced and a small change in length of two target probes. |
| A3 | Two reference probes are replaced and the length of several probes adjusted. |
| A2 | QDX2 fragments have been added. |
| A1 | First release. |

Implemented changes in the product description

Version A4-02 — 05 July 2022 (02P)

- Corrected the exon numbering for the *KCNQ1* probes according to LRG_287 and added old exon numbering between brackets to Table 2e.
- Various minor textual or layout changes.

Version A4-01 — 24 March 2020 (02P)

- Product description rewritten and adapted to a new template.
- Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2).
- Ligation sites of the probes targeting the *CHRNA2*, *NHLRC1*, *EPM2A*, *KCNQ3*, and *KCNQ1* genes updated according to new versions of the NM_ reference sequences.
- Warning removed in Tables for 214 nt probe 06609-L06167, 240 nt probe 06610-L10371 and 400 nt probe 06608-L06166.
- Exon numbering has been adjusted for 274 nt *KCNQ3* probe and 400 nt *CHRNA4* probe from exon 1 to Upstream.

Version 13 – 15 September 2017 (55)

- Warning added in Table 1, 160 nt probe 06617-L29224, 214 nt probe 06609-L06167, 238 nt probe 06610-L29157 274 nt probe 06593-L06151, 293 nt probe 06615-L06173, 400 nt probe 06608-L06166, and 454 nt probe 08195-L08089.

Version 12 – 25 May 2016 (55)

- Product description adapted to a new product version (version number changed, lot number added, small changes in Table 1 and Table 2, new picture included).
- References on page 2 updated.
- Exon numbering of the *KCNQ3* and *KCNQ1* genes has been changed in Table 1 and Table 2.

Version 11 – 15 July 2015 (54)

- Figure based on the use of old MLPA buffer (replaced in December 2012) removed.

Version 10 (48)

- Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added.

More information: www.mlpa.com; www.mlpa.eu



MRC-Holland bv; Willem Schoutenstraat 1
1057 DL, Amsterdam, The Netherlands

| | |
|--------|---|
| E-mail | info@mlpa.com (information & technical questions); order@mlpa.com (orders) |
| Phone | +31 888 657 200 |