

Product Description

SALSA® MLPA® Probemix P489-A2 BARD1

To be used with the MLPA General Protocol.

Version A2

As compared to version A1, one reference probe has been replaced. For complete product history see page 7.

Catalogue numbers

- **P489-025R:** SALSA® MLPA® Probemix P489 BARD1, 25 reactions
- **P489-050R:** SALSA® MLPA® Probemix P489 BARD1, 50 reactions
- **P489-100R:** SALSA® MLPA® Probemix P489 BARD1, 100 reactions

SALSA® MLPA® Probemix P489 BARD1 (hereafter: P489 BARD1) is to be used in combination with:

1. SALSA® MLPA® Reagent Kit (Cat. No: EK1-FAM, EK1-CY5, EK5-FAM, EK5-CY5, EK20-FAM),
2. Data analysis software Coffalyser.Net™ (Cat. No: n.a.)

Volumes and ingredients

Volumes			Ingredients
P489-025R	P489-050R	P489-100R	
40 µl	80 µl	160 µl	Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA

The MLPA probemix is not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. A Safety Data Sheet (SDS) is not required for this product: none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).

Storage and handling

Recommended storage conditions		
--------------------------------	---	---

A shelf life of until the expiry date is guaranteed, when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. This product should not be exposed to more than 25 freeze-thaw cycles. Do not use the product if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations.

Certificate of Analysis

Information regarding quality tests and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

SALSA® MLPA® Probemix P489 BARD1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *BARD1* gene, which is associated with hereditary breast and ovarian cancer, as well as neuroblastoma.

The *BARD1* (BRCA1 associated RING domain 1) gene is transcribed from chromosome 2q35 and functions as tumour suppressor by playing an essential role in BRCA1-mediated double-strand DNA repair. *BARD1* shares

homology with the two most conserved regions of the cancer susceptibility gene *BRCA1*. Both *BARD1* and *BRCA1* harbour two BRCA1 C-terminal repeats (BRCT) and an N-terminal RING motif (Tarsounas et al. 2020). Heterodimerization of the RING motives of both proteins is essential for the E3 ubiquitin ligase activity of the BRCA1-BARD1 complex. Moreover, the BRCA1-BARD1 heterodimer facilitates DNA repair via the homologous recombination pathway by interacting with DNA repair proteins such as BRCA2 and RAD51.

Autosomal dominant mutations in the genes *BRCA1* and *BRCA2* are the most frequent cause for familial breast cancer. Other breast cancer susceptibility genes include *PALB2*, *CHEK2*, *ATM*, and *BARD1*. While most pathogenic variants of *BARD1* are caused by small point mutations, copy number variations have also been reported (Sabatier et al. 2010). Moreover, copy number variations of *BARD1* have been related to congenital conditions such as hypospadias and congenital heart defects as well as developmental delay (Cimmino et al. 2017).

This product is not CE/FDA registered for use in diagnostic procedures. The SALSA® MLPA® technique is covered by US patent 6,955,901 and corresponding patents outside the US. The purchase of this product includes a license to use only this amount of product solely for the purchaser's own use.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <https://www.ncbi.nlm.nih.gov/gene>

For NM_ mRNA reference sequences: <https://www.ncbi.nlm.nih.gov/nucore?db=nucleotide>

Matched Annotation from NCBI and EMBL-EBI (MANE): <https://www.ncbi.nlm.nih.gov/refseq/MANE>

Tark – Transcript Archive: <https://tark.ensembl.org>

Exon numbering

The *BARD1* exon numbering used in this P489-A2 BARD1 product description is the exon numbering from the LRG_297 sequence. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

P489-A2 BARD1 contains 25 MLPA probes with amplification products between 138 and 346 nucleotides (nt). This includes 13 probes for the *BARD1* gene, including at least one probe for every exon. Moreover, two probes upstream of the *BARD1* gene are included. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of MLPA (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation using 16 different DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample type or the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, free from impurities known to affect MLPA reactions. MRC Holland has tested and can recommend the following extraction methods:

- QIAGEN Autopure LS (automated) and QIAamp DNA mini/midi/maxi kit (manual)
- Promega Wizard Genomic DNA Purification Kit (manual)
- Salting out (manual)

All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of different reference samples from unrelated individuals should be included in each MLPA experiment for data normalisation. Reference samples should be derived from individuals who are from families without a history of hereditary cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore deviations to the indicated copy number alteration (CNA) findings might occur.

Data analysis

Coffalyser.Net should be used for data analysis in combination with the appropriate lot-specific Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	FR = 0
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication/gain	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication/gain	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of Coffalyser.Net (Calculations, cut-offs and interpretation remain unchanged.) Please note that Coffalyser.Net also shows arbitrary borders as part of the statistical analysis of

results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. single nucleotide variants (SNVs), point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination in the DNA sample) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *BARD1* gene are small (point) mutations, none of which will be detected by using P489 BARD1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism in sequence data indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

BARD1 mutation database

<https://databases.lovd.nl/shared/genes/BARD1>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <https://varnomen.hgvs.org>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *BARD1* exons 5 and 7 but not exon 6) to MRC Holland: info@mrcholland.com.

Table 1. P489-A2 BARD1

Length (nt)	MLPA probe	Chromosomal position (hg18) ^a	
		Reference	<i>BARD1</i>
64-105	Control fragments – see table in probemix content section for more information		
138	Reference probe 16316-L29209	3q	
147	BARD1 probe 22759-L32103		Exon 6
153 –	BARD1 probe 22760-L32104		Upstream
165	Reference probe 09632-L09917	17q	
173	BARD1 probe 22761-L32105		Exon 10
178	BARD1 probe 22762-L32106		Exon 7
190	Reference probe 11368-L12093	4q	
198	BARD1 probe 22763-L32107		Exon 11
206	BARD1 probe 22764-L32180		Exon 3
214	Reference probe 10730-L30523	6p	
220	BARD1 probe 22765-L32109		Exon 1
228	BARD1 probe 22766-L32110		Exon 9
239	Reference probe 15025-L16774	11p	
246	BARD1 probe 22767-L32111		Exon 4
256 –	BARD1 probe 22768-L32112		Upstream
265	Reference probe 11438-L12168	1q	
274	Reference probe 08053-L07834	5p	
283	BARD1 probe 22769-L32113		Exon 11
292	BARD1 probe 22770-L32114		Exon 5
301	Reference probe 16436-L26274	18q	
310	Reference probe 18380-L25673	10q	
319	BARD1 probe 22771-L32115		Exon 2
328	BARD1 probe 22772-L32116		Exon 4
337	BARD1 probe 22773-L32117		Exon 8
346 *	Reference probe 20766-L28668	1q	

^a See section Exon numbering on page 2 for more information.

* New in version A2.

– Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

The probe lengths in the table above may vary slightly depending on the capillary electrophoresis machine settings. Please see the most up-to-date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. *BARD1* probes arranged according to chromosomal location

Length (nt)	MLPA probe	<i>BARD1</i> exon ^a	Ligation site ^b NM_000465.4	Partial sequence ^c (24 nt adjacent to ligation site)	Distance to next probe
256 -	22768-L32112	Upstream	1685 nt before exon 1, reverse	GTTGAATCATTG-AACTACCTACCC	0.4 kb
153 -	22760-L32104	Upstream	1269 nt before exon 1	AAAGAGCAGAGA-AAACTTTAAGG	1.5 kb
		<i>start codon</i>	115-117 (Exon 1)		
220	22765-L32109	Exon 1	193-194	CGCCCCCATGG-AACCGGATGGTC	12.4 kb
319	22771-L32115	Exon 2	311-312	TTTAGGAGGATG-TGAGCACATCTT	4.7 kb
206	22764-L32180	Exon 3	436-437	TGATTCAACTTT-GTAGTAAGCTTC	10.9 kb
246	22767-L32111	Exon 4	515-516	AAGTTTGTAA-TGATGCAGGAAA	0.6 kb
328	22772-L32116	Exon 4	1093-1094	GCCATCACAATA-GACTTTCCAGTC	11.6 kb
292	22770-L32114	Exon 5	1476-1477	GGAAGTGATCCA-AATGTTAAAGAC	1.7 kb
147	22759-L32103	Exon 6	1575-1574, reverse	CCGGTGGTGTTC-ACCAATGCCTTA	15.1 kb
178	22762-L32106	Exon 7	1740-1739, reverse	GGTAGCAGCAAT-AGCGATTTCATA	6.7 kb
337	22773-L32117	Exon 8	1861-1862	CAGAACAACAGA-AAATGCTCAGTG	0.7 kb
228	22766-L32110	Exon 9	1967-1968	TCAAAGTACCTT-GAAGTGATGCT	14.6 kb
173	22761-L32105	Exon 10	2055-2054, reverse	TACTTTCTTCC-TGTTCCACATACT	1.6 kb
198	22763-L32107	Exon 11	2258-2257, reverse	TATTGATGGTCT-GAGTCACGTCAC	1.7 kb
283	22769-L32113	Exon 11	3939-3940	ATATATGAACAA-AATAATAACCTT	
		<i>stop codon</i>	2446-2448 (Exon 11)		

^a See section Exon numbering on page 2 for more information.

^b Ligation sites are relative to the start of the NM_ sequence, and not relative to the coding sequence.

^c Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

- Flanking probe. Included to help determine the extent of a deletion/duplication. Copy number alterations of only the flanking or reference probes are unlikely to be related to the condition tested.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Related products

P002 BRCA1	Contains probes for the <i>BRCA1</i> gene, involved in breast and ovarian cancer.
P087 BRCA1 Confirmation	Contains probes for the <i>BRCA1</i> gene. It can be used to confirm the results obtained with P002.
P045 BRCA2/CHEK2	Contains probes for the <i>BRCA2</i> and <i>CHEK2</i> genes, involved in breast and ovarian cancer.
P090 BRCA2	Contains the same probes for the <i>BRCA2</i> gene as P045.
P077 BRCA2 Confirmation	Contains probes for the <i>BRCA2</i> gene. It can be used to confirm the results obtained with P045 or P090.
P041/P042 ATM	Contain probes for the <i>ATM</i> gene, involved in breast cancer and Ataxia Telangiectasia.
P056 TP53	Contains probes for <i>TP53</i> , involved in Li-Fraumeni syndrome.
P260 PALB2-RAD50-RAD51C-RAD51D	Contains probes for the <i>PALB2</i> , <i>RAD50</i> , <i>RAD51C</i> , and <i>RAD51D</i> genes, involved in breast and ovarian cancer.

References

- Cimmino F et al. (2017). Dualistic role of *BARD1* in cancer. *Genes (Basel)*. 8:e375.
- Sabatier R et al. (2010). *BARD1* homozygous deletion, a possible alternative to *BRCA1* mutation in basal breast cancer. *Gene Chromosomes Canc.* 49:1143-51.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.

- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Tarsounas M et al. (2020). The antitumorigenic roles of BRCA1-BARD1 in DNA repair and replication. *Nat Rev Moll Cell Biol.* 21:289-99.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using P489 BARD1

- Benito-Sánchez B et al. (2022). Apparent regional differences in the spectrum of BARD1 pathogenic variants in Spanish population and importance of copy number variants. *Sci Rep*, 12(1), 8547.
- Carrera S et al. (2023). Germline heterozygous exons 8–11 pathogenic BARD1 gene deletion reported for the first time in a family with suspicion of a hereditary colorectal cancer syndrome: more than an incidental finding? *Hered Cancer Clin Pract*, 21(1), 2.

P489 product history	
Version	Modification
A2	One reference probe has been replaced.
A1	First release.

Implemented changes in the product description
Version A2-01 – 23 August 2024 (04P) <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2). - New references added to section “Selected publications using P489 BARD1”.
Version A1-02 – 19 January 2021 (04P) <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template.
Version A1-01 – 05 June 2020 (02P) <ul style="list-style-type: none"> - Not applicable, new document.

More information: www.mrcholland.com ; www.mrcholland.eu	
	MRC Holland BV; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200