

Product Description

SALSA® MLPA® Probemix P433-A2 ARID1A-ARID1B

To be used with the MLPA General Protocol.

Version A2

As compared to version A1, several probes have a change in length but no change in the sequence targeted. For complete product history see page **Error! Bookmark not defined.**

Catalogue numbers

- **P433-025R:** SALSA® MLPA® Probemix P433 ARID1A-ARID1B, 25 reactions
- **P433-025R:** SALSA® MLPA® Probemix P433 ARID1A-ARID1B, 50 reactions
- **P433-100R:** SALSA® MLPA® Probemix P433 ARID1A-ARID1B, 100 reactions

SALSA® MLPA® Probemix P433 ARID1A-ARID1B (hereafter: P433 ARID1A-ARID1B) is to be used in combination with:

1. SALSA® MLPA® Reagent Kit (Cat. No: EK1-FAM, EK1-CY5, EK5-FAM, EK5-CY5, EK20-FAM),
2. Data analysis software Coffalyser.Net™ (Cat. No: n.a.)

Volumes and ingredients

Volumes			Ingredients
P433-025R	P433-050R	P433-100R	
40 µl	80 µl	160 µl	Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA

The MLPA probemix is not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. A Safety Data Sheet (SDS) is not required for this product: none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).

Storage and handling

Recommended storage conditions		
--------------------------------	---	---

A shelf life of until the expiry date is guaranteed, when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. This product should not be exposed to more than 25 freeze-thaw cycles. Do not use the product if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations.

Certificate of Analysis

Information regarding quality tests and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

SALSA® MLPA® Probemix P433 ARID1A-ARID1B is a **research use only (RUO)** assay for the detection of deletions or duplications in the *ARID1A* and *ARID1B* genes, which are associated with various tumour types, Coffin-Siris syndrome (OMIM #614607) and non-syndromic intellectual disability (OMIM #614562).

ARID1A (AT-rich interactive domain 1A) and ARID1B (AT-rich interactive domain 1B) proteins share 80% amino acid homology and are evolutionarily conserved. These proteins are integral components of the SWI/SNF chromatin modelling complex, which mediates epigenetic regulation in critical cellular processes such as cell differentiation, proliferation and DNA repair (Madan V et al. 2023).

Alterations in *ARID1A* and *ARID1B* genes are typically loss-of-function changes and they can exhibit bi-allelic inactivation or loss of protein expression, consistent with a tumour suppressor mechanism (Garraway and Lander, 2013). *ARID1A* is frequently mutated across a wide variety of human cancers, whereas alterations in *ARID1B* have been detected at lower frequencies. *ARID1A* and *ARID1B* genetic alterations have been identified in 11% of childhood neuroblastoma and are associated with early treatment failure and decreased survival (Sausen et al. 2012). Furthermore, *ARID1A* inactivating mutations have been detected in approximately 50% of ovarian clear cell carcinomas (Jones et al. 2010, Wiegand et al. 2010, Caumanns JJ et al. 2018), in up to 17% of hepatocellular carcinomas (Guichard et al. 2012, Fujimoto et al. 2012, Huang et al. 2012), in 10-29% of gastric cancer cases (Wang et al. 2011, Zang et al. 2012) and in 37% of breast cancers (Cornen et al. 2012). Importantly, it has been shown that ARID1A-deficiency sensitizes cancer cells to PARP inhibitor therapies, providing new treatment possibilities for ARID1A-mutant tumours (Shen et al. 2015; Park et al. 2019; Yu et al. 2023).

Genetic alterations in *ARID1A* (<5%) and *ARID1B* (~37%) have also been detected in Coffin-Siris syndrome and non-syndromic intellectual disability cases. *ARID1A* microduplications cause intellectual disability with recognizable syndromic features (Bidart et al. 2017), whereas intragenic duplications and microdeletions of *ARID1B* have been described in patients with severe intellectual disability (Hoyer et al. 2012). *ARID1B* alterations have also been found in patients with Coffin-Siris syndrome (Santen et al. 2012, Tsurusaki et al. 2012) and in patients affected with either syndromic or non-syndromic short stature (Yu et al. 2015).

This product is not CE/FDA registered for use in diagnostic procedures. The SALSA® MLPA® technique is covered by US patent 6,955,901 and corresponding patents outside the US. The purchase of this product includes a license to use only this amount of product solely for the purchaser's own use.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <https://www.ncbi.nlm.nih.gov/gene>

For NM_ mRNA reference sequences: <https://www.ncbi.nlm.nih.gov/nuccore?db=nucleotide>

Matched Annotation from NCBI and EMBL-EBI (MANE): <https://www.ncbi.nlm.nih.gov/refseq/MANE>

Exon numbering

The *ARID1A* and *ARID1B* exon numbering used in this P433-A2 ARID1A-ARID1B product description is the exon numbering from the MANE project (release version 1.3), based on MANE Select transcripts NM_006015.6 and NM_001374828.1, respectively. The *ARID1B* exon numbering has changed; the exon numbering used in previous versions of this product description can be found in between brackets in Table 2. From description version A2-02 onwards, we have adopted the NCBI exon numbering that is present in the MANE Select NM_ sequence for this gene. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

P433-A2 ARID2A-ARID2B contains 56 MLPA probes with amplification products between 127 and 504 nucleotides (nt). This includes 22 probes for the *ARID2A* and 22 probes for *ARID1B* gene. In addition, 12 reference probes are included that detect autosomal chromosomal locations that target relatively stable copy number regions in various cancer types. Complete probe sequences and the identity of the genes detected by the reference probes are available in Table 3 and online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below).

More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of MLPA (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com). More information on the use of MLPA in tumour applications can be found in Hömig-Hölzel and Savola (2012).

MLPA technique validation

Internal validation using 16 different DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample type or the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA, which includes DNA derived from formalin-fixed, paraffin-embedded (FFPE) tissues, free from impurities known to affect MLPA reactions. MRC Holland has tested and can recommend the following extraction methods:

- QIAGEN Autopure LS (automated) and QIAamp DNA mini/midi/maxi kit (manual)
- Promega Wizard Genomic DNA Purification Kit (manual)
- Salting out (manual)

All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol. More information on the use of FFPE tissue samples for MLPA can be found in Atanesyan et al. (2017).

Reference samples

A sufficient number (≥ 3) of different reference samples from healthy individuals without history of cancer should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples for germline analysis should be derived from different unrelated individuals who are from families without a history of Coffin-Siris syndrome and intellectual disability, and for tumour analysis from healthy individuals without history of cancer. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ (<https://www.dsmz.de>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. Samples mentioned below from the Coriell Institute have been tested with P433-A2 ARID1A-ARD1B probemix at MRC Holland and can be used as positive control samples. The quality of cell lines can change; therefore deviations to the indicated copy number alteration (CNA) findings might occur.

Sample name	Source	Chromosomal position of CNA*	Altered target genes in P433-A2 ARID1A-ARID1B	Expected copy number alteration
NA07994	Coriell Institute	6q25.3	ARID1B	Heterozygous duplication
NA06802	Coriell Institute	6q25.3	ARID1B	Heterozygous deletion

* Indicated chromosomal bands accommodate genes targeted by MLPA probes, however, the whole extent of copy number alteration (CNA) present in this cell line cannot be determined by this P433-A2 ARID1A-ARID1B Probemix.

Data analysis

Coffalyser.Net should be used for data analysis in combination with the appropriate lot-specific Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 . When this criterion is fulfilled, the following cut-off values for the final ratio (FR) of the probes can be used to interpret MLPA results for autosomal chromosomes or pseudo-autosomal regions:

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	$FR = 0$
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication/gain	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication/gain	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of Coffalyser.Net (Calculations, cut-offs and interpretation remain unchanged.) Please note that Coffalyser.Net also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

Please note that these above mentioned final ratios are only valid for germline testing. Final ratios are affected both by percentage of tumour cells and by possible subclonality.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic or subclonal cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. single nucleotide variants (SNVs), point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination in the DNA sample) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *ARID1A* or *ARID1B* genes. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.

- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of germline genetic defects in the *ARID1A* and *ARID1B* genes are small (point) mutations, none of which will be detected by using P433 ARID1A-ARID1B. Likewise, in many tumour samples genetic alterations in the *ARID1A* and *ARID1B* genes are small (point) mutations, none of which will be detected by using P433 ARID1A-ARID1B.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.
- MLPA analysis on tumour samples provides information on the *average* situation in the cells from which the DNA sample was purified. Gains or losses of genomic regions or genes may not be detected if the percentage of tumour cells is low. In addition, subclonality of the aberration affects the final ratio of the corresponding probe. Furthermore, there is always a possibility that one or more reference probes *do* show a copy number alteration in a patient sample, especially in solid tumours with more chaotic karyotypes.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism in sequence data indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

ARID1A and ARID1B mutation database

We strongly encourage users to deposit positive results in the COSMIC (<http://cancer.sanger.ac.uk/cosmic>) and LOVD (<https://databases.lovd.nl>) mutation databases. Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report false positive results due to SNVs and unusual results (e.g., a duplication of *ARID1A* exons 4 and 6 but not exon 5) to MRC Holland: info@mrcholland.com.

Table 1. P433-A2 ARID1A-ARID1B

Length (nt)	MLPA probe	Chromosomal position (hg18) ^a		
		Reference	ARID1A, 1p36.11	ARID1B, 6q25.3
64-105	Control fragments – see table in probemix content section for more information			
127	Reference probe 18946-L24593	5q31		
134 « Δ	ARID1A probe 18899-L24494		Exon 1	
142	ARID1A probe 18900-L24495		Exon 7	
148	ARID1B probe 19238-L25319			Exon 13
154	ARID1B probe 18902-L24497			Exon 20
160 « ¥ Δ	ARID1B probe 18903-L32281			Exon 1
166	Reference probe 12741-L21124	21q22		
172 «	ARID1A probe 18904-L25189		Exon 1	
178	ARID1A probe 18905-L24500		Exon 5	
184	ARID1B probe 18906-L24501			Exon 15
190	ARID1B probe 18907-L24502			Exon 18
196	ARID1B probe 18908-L24503			Exon 9
202 «	ARID1B probe 18909-L24504			Exon 1
208	Reference probe 11546-L13329	19q13		
214 ¥ Δ	ARID1A probe 18910-L32453		Exon 8	
220	ARID1A probe 18911-L24506		Exon 2	
226	ARID1A probe 18912-L24507		Exon 10	
232 «	ARID1B probe 18913-L25477			Exon 8
238	Reference probe 08070-L07851	9p13		
244	ARID1B probe 18914-L24509			Exon 5
250	ARID1A probe 18915-L24510		Exon 4	
256	ARID1B probe 18916-L24511			Exon 10
262	ARID1B probe 18917-L24512			Exon 6
268 Ø	ARID1B probe 18918-L24513			Intron 2
274	Reference probe 17873-L22132	2p21		
281	ARID1B probe 18919-L25190			Exon 12
287	ARID1A probe 18920-L25191		Exon 16	
292	ARID1A probe 18921-L25192		Exon 3	
300	ARID1B probe 18922-L24517			Exon 7
308	ARID1A probe 18923-L24518		Exon 15	
316	ARID1A probe 18924-L25193		Exon 19	
322	ARID1A probe 18925-L25194		Exon 9	
328	ARID1A probe 18926-L25195		Exon 17	
337	Reference probe 03959-L20754	11p11		
346	ARID1A probe 18927-L24522		Exon 20	
353	ARID1A probe 18928-L24523		Exon 12	
359	ARID1B probe 18929-L25196			Exon 2
367	Reference probe 15524-L25197	16q13		
373	ARID1B probe 18930-L24525			Exon 14
380	ARID1B probe 18931-L24526			Exon 4
391	ARID1A probe 18932-L24527		Exon 14	
398	ARID1B probe 18933-L25198			Exon 19
404 ¥	Reference probe 21393-L25199	9q21		
409	ARID1A probe 18934-L25200		Exon 18	
417	ARID1B probe 18935-L25201			Exon 17
423	ARID1B probe 18936-L25202			Exon 3
427	ARID1A probe 18937-L25203		Exon 11	
436	ARID1A probe 18938-L25229		Exon 13	
445	Reference probe 13810-L15304	5q14		
451 Ø	ARID1B probe 18939-L24534			Intron 4
460	ARID1A probe 18940-L25206		Exon 18	

Length (nt)	MLPA probe	Chromosomal position (hg18) ^a		
		Reference	ARID1A, 1p36.11	ARID1B, 6q25.3
466	Reference probe 00979-L25207	10p14		
472	ARID1A probe 18941-L24536		Exon 6	
481	ARID1B probe 18942-L24537			Exon 16
490	Reference probe 14909-L17529	18p11		
504	Reference probe 09870-L19465	2p15		

^a See section Exon numbering on page 2 for more information.

¥ Changed in version A2. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Ø Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

The probe lengths in the table above may vary slightly depending on the capillary electrophoresis machine settings. Please see the most up-to-date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. P433-A2 target probes arranged according to chromosomal location

Table 2a. ARID1A, 1p36.11

Length (nt)	MLPA probe	Gene exon ^a	Ligation site ^b NM_006015.6	Partial sequence ^c (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	390-392 (Exon 1)		
134 « Δ	18899-L24494	Exon 1	384-385	ACGAGACAGCGG-GGATCATGGCCG	0.6 kb
172 «	18904-L25189	Exon 1	1024-1025	CAACCACCAGTA-CAACTCCTACTA	32.6 kb
220	18911-L24506	Exon 2	1537-1538	GCCATCCAGTCC-AATGGATCAGAT	1.9 kb
292	18921-L25192	Exon 3	2189-2190	TTCCCTCCACCG-CAGGTAAGATAT	1.2 kb
250	18915-L24510	Exon 4	2286-2287	ACCTGAGCCTTC-AGTCAAGACCCT	28.2 kb
178	18905-L24500	Exon 5	2385-2386	GCACATCAGGGA-TTCCAGCAGCC	0.5 kb
472	18941-L24536	Exon 6	2638-2639	TGCCCAAGATCG-AGGTGAGAGCCT	0.9 kb
142	18900-L24495	Exon 7	25 nt after Exon 7	CCAGTTAGGAGT-AGATACGGGTGA	0.7 kb
214 ¥ Δ	18910-L32453	Exon 8	2867-2866 reverse	CCAGCCATGCCT-GCACTGGGGTAG	3.3 kb
322	18925-L25194	Exon 9	3182-3183	CCTCCTTATGGA-CAAGGGATTAAT	0.3 kb
226	18912-L24507	Exon 10	3370-3371	GACAGAATCCAA-ATCCAAGGTAGT	1.3 kb
427	18937-L25203	Exon 11	3403-3404	TACTACAACCAA-TGAGAAGATCAC	3.4 kb
353	18928-L24523	Exon 12	3656-3657	AGCAGTGCTGCC-AGCTCCTTGA	1.4 kb
436	18938-L25229	Exon 13	3924-3923 reverse	GCCTTACCTCAT-GCCTGGCAATGG	0.4 kb
391	18932-L24527	Exon 14	6 nt after Exon 14	GGAAAGGTGACT-GATCTGATTGCT	0.3 kb
308	18923-L24518	Exon 15	16 nt before Exon 15	TGTTACCGCTT-GCCTTTCTACGC	0.3 kb
287	18920-L25191	Exon 16	4284-4283 reverse	AGTGCTCATGTT-TCCCTCAGGCC	0.3 kb
328	18926-L25195	Exon 17	4485-4486	ATCAACAGCAAC-AGCAGGTGAGGA	0.6 kb
409	18934-L25200	Exon 18	4691-4692	TATGGCAATGCC-TATCCTGCCACT	0.5 kb
460	18940-L25206	Exon 18	5204-5205	AGCAAGTCTCCA-TTCCTGCACTCT	0.6 kb
316	18924-L25193	Exon 19	5453-5454	ACATGGGCATTA-GATACCATCAAC	3.4 kb
346	18927-L24522	Exon 20	5546-5547	CTTGTAAGATAT-TTCCGACGATGC	-
		<i>stop codon</i>	7245-7247 (Exon 20)		

Table 2b. *ARID1B*, 6q25.3

Length (nt)	MLPA probe	Gene exon ^a	Ligation site ^b NM_001374828.1	Partial sequence ^c (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	304-306 (Exon 1)		
202 «	18909-L24504	Exon 1	1323-1324	CAACATGGCGGA-CAACAAAGCCCC	0.5 kb
160 « ¥ Δ	18903-L32281	Exon 1	1848-1849	ATGATGCGGAGC-TACGGCGGCAGC	50.1 kb
359	18929-L25196	Exon 2	2193-2194	CCGTACCCAGGA-GGTTCTATGGC	42.3 kb
268 Ø	18918-L24513	Intron 2 (Exon 3)	29.7 kb before Exon 3 (NM_001374820.1; 2321-2322)	GAAAGAAACATT-CTGGTTGGCAAG	29.9 kb
423	18936-L25202	Exon 3 (4)	2433-2432 reverse	CTCACCTGCTGC-GGCTGGTACCTC	34.0 kb
380	18931-L24526	Exon 4 (5)	2545-2544 reverse	TCTTACTGGTAA-ACTTGATGGTCT	101.3 kb
451 Ø	18939-L24534	Intron 4 (5)	47.8 kb before Exon 5 (NM_001363725.2; 128-129)	GATGCAGCCTTT-GGACTCAAGGTG	48.0 kb
244	18914-L24509	Exon 5 (6)	2743-2744	TTGGCTCTCCTG-TAGGAAGCAACC	25.7 kb
262	18917-L24512	Exon 6 (7)	2834-2835	GAGCCAGTCAGA-ATCCAGTTCCCA	22.7 kb
300	18922-L24517	Exon 7 (8)	21 nt after Exon 7	AAGCTTCTCCAA-AATGCATGGCAG	15.7 kb
232 «	18913-L25477	Exon 8 (9)	3330-3331	CCAAGTGTGAAC-CGTAAAGGCACAG	18.2 kb
196	18908-L24503	Exon 9 (10)	3452-3453	CTCTCCCTACAG-CCAGCCCATGAA	7.0 kb
256	18916-L24511	Exon 10 (11)	3582-3583	GGTGAATCCAAA-CTGCCCTGCCT	7.0 kb
281	18919-L25190	Exon 12 (13)	3879-3880	GAGCCAGAGAGA-AAGCTCTGGGTC	3.4 kb
148	19238-L25319	Exon 13 (14)	4219-4220	AGCCGCCATCTC-CTGGTAAGTGGC	5.3 kb
373	18930-L24525	Exon 14 (15)	4290-4291	ATGGCAGAGGTT-CCAGGTGACCTG	0.4 kb
184	18906-L24501	Exon 15 (16)	4410-4411	GTGAGTGATTCA-TCCTTCCCGAAA	6.2 kb
481	18942-L24537	Exon 16 (17)	4608-4609	ATGTACAACCAA-AGTCCCTCCGGA	2.6 kb
417	18935-L25201	Exon 17 (18)	4709-4710	TGGGCAGCAGTA-TCCAGGCCAAGG	2.2 kb
190	18907-L24502	Exon 18 (19)	5072-5073	TATGTGGGCAGC-ACGCAATGATAT	3.0 kb
398	18933-L25198	Exon 19 (20)	5672-5673	TGATGACAGCAC-TGTTGCTACTTT	3.6 kb
154	18902-L24497	Exon 20 (21)	7104-7105	GGGGACGCACTA-GCAGCAAGGGCC	-
		<i>stop codon</i>	7420-7422 (Exon 20)		

^a See section Exon numbering on page 2 for more information.

^b Ligation sites are relative to the start of the NM_ sequence, and not relative to the coding sequence.

^c Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

¥ Changed in version A2. Minor alteration, no change in sequence detected.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Δ More variable. This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution.

Ø Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridisation and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 3. Reference probes arranged according to chromosomal location

Length (nt)	MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence (24 nt adjacent to ligation site)	Location (hg18) in kb
274	17873-L22132	<i>PPM1B</i>	2p21	AGCAGAAAATCA-TTAGCATTTCCC	02-044,313
504	09870-L19465	<i>PEX13</i>	2p15	TGAGGATGACCA-TGTAGTTGCCAG	02-061,126
445	13810-L15304	<i>ADGRV1</i>	5q14	AAAGGGAGTTTA-GCCTGGCATTAG	05-090,060
127	18946-L24593	<i>IL4</i>	5q31	ATCGACACCTAT-TAATGGGTCTCA	05-132,038

Length (nt)	MLPA probe	Gene	Chromosomal band (hg18)	Partial sequence (24 nt adjacent to ligation site)	Location (hg18) in kb
238	08070-L07851	<i>DNAI1</i>	9p13	AGGGATCTGTTT-GCAGTGGGATAT	09-034,491
404	21393-L25199	<i>PCSK5</i>	9q21	GATGAGCTGGAA-TATGATGACGAG	09-078,164
466	00979-L25207	<i>UPF2</i>	10p14	TGCCATTCTTTT-GCATCTCAAAAG	10-012,019
337	03959-L20754	<i>MYBPC3</i>	11p11	AGGCTGACTACA-GCTTTGTGCCCCG	11-047,319
367	15524-L25197	<i>SLC12A3</i>	16q13	CACCCGGAACCT-CAGCCTGATGAT	16-055,478
490	14909-L17529	<i>RNMT</i>	18p11	TACAATGAACCT-CAGGAAGTTGGT	18-013,724
208	11546-L13329	<i>PRPF31</i>	19q13	TGAGCTTACTGA-TCATGATAGGAC	19-059,314
166	12741-L21124	<i>RIPK4</i>	21q22	AAGCCAAGAAGA-TGGAGATGGCCA	21-042,050

Complete probe sequences are available at www.mrcholland.com.

References

- Atanesyan L et al. (2017). Optimal fixation conditions and DNA extraction methods for MLPA analysis on FFPE tissue-derived DNA. *Am J Clin Pathol.* 147:60-8.
- Balbás-Martínez C et al. (2013). ARID1A alterations are associated with FGFR3-wild type, poor-prognosis, urothelial bladder tumors. *PLoS ONE.* 8:e62483.
- Bidart M et al. (2017). Microduplication of the ARID1A gene causes intellectual disability with recognizable syndromic features. *Genet Med.* 19:701-10.
- Caumanns JJ et al. (2018). ARID1A mutant ovarian clear cell carcinoma: A clear target for synthetic lethal strategies. *Biochim Biophys Acta Rev Cancer.* 1870:176-84.
- Cornen S et al. (2012). Mutations and deletions of ARID1A in breast tumors. *Oncogene.* 31:4255-6.
- Fujimoto A et al. (2012). Whole-genome sequencing of liver cancers identifies etiological influences on mutation patterns and recurrent mutations in chromatin regulators. *Nat Genet.* 44:760-4.
- Garraway LA and Lander S (2013). Lessons from the cancer genome. *Cell.* 153:17-37.
- Guichard C et al. (2012). Integrated analysis of somatic mutations and focal copy-number changes identifies key genes and pathways in hepatocellular carcinoma. *Nat Genet.* 44:694-8.
- Hömig-Hölzel C and Savola S. (2012). Multiplex ligation-dependent probe amplification (MLPA) in tumor diagnostics and prognostics. *Diagn Mol Pathol.* 21:189-206.
- Hoyer J et al. (2012). Haploinsufficiency of ARID1B, a member of the SWI/SNF-a chromatin-remodeling complex, is a frequent cause of intellectual disability. *Am J Hum Genet.* 90:565-72.
- Huang J et al. (2012). Exome sequencing of hepatitis B virus-associated hepatocellular carcinoma. *Nat Genet.* 44:1117-21.
- Jones S et al. (2010). Frequent mutations of chromatin remodeling gene ARID1A in ovarian clear cell carcinoma. *Science.* 330:228-31.
- Madan V et al. (2023). Dissecting the role of SWI/SNF component ARID1B in steady-state hematopoiesis. *Blood Adv.* 7:6553-66.
- Park Y et al. (2019). Loss of ARID1A in Tumor Cells Renders Selective Vulnerability to Combined Ionizing Radiation and PARP Inhibitor Therapy. *Clin Cancer Res.* 25:5584-94.
- Sausen M et al. (2012). Integrated genomic analyses identify ARID1A and ARID1B alterations in the childhood cancer neuroblastoma. *Nat Genet.* 45:12-7.
- Santen GWE et al. (2012). Mutations in SWI/SNF chromatin remodelling complex gene ARID1B cause Coffin-Siris syndrome. *Nat Genet.* 44:379-80.
- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.
- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Shen J et al. (2015). ARID1A Deficiency Impairs the DNA Damage Checkpoint and Sensitizes Cells to PARP Inhibitors. *Cancer Discov.* 5:752-67.
- Tsurusaki Y et al. (2012). Mutations affecting components of the SWI/SNF complex cause Coffin-Siris syndrome. *Nat Genet.* 44:376-8.

- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.
- Wang K et al. (2011). Exome sequencing identifies frequent mutation of ARID1A in molecular subtypes of gastric cancer. *Nat Genet.* 43:1219-23.
- Wiegand KC et al. (2010). ARID1A mutations in endometriosis-associated ovarian carcinomas. *N Engl J Med.* 363:1532-43.
- Yu Y. et al. (2015). De novo mutations in ARID1B associated with both syndromic and non-syndromic short stature. *BMC Genomics.* 16:701.
- Yu ZC et al. (2023). Temozolomide Sensitizes ARID1A-Mutated Cancers to PARP Inhibitors. *Cancer Res.* 83:2750-2762.
- Zang ZJ et al. (2012). Exome sequencing of gastric adenocarcinoma identifies recurrent somatic mutations in cell adhesion and chromatin remodeling genes. *Nat Genet.* 44:570-4.

Selected publications using P433 ARID1A-ARID1B

- Yu Y. et al. (2015). De novo mutations in ARID1B associated with both syndromic and non-syndromic short stature. *BMC Genomics.* 16:701.

P433 product history	
Version	Modification
A2	Several probes have a change in length but no change in the sequence targeted.
A1	First release.

Implemented changes in the product description
<p>Version A2-02 – 24 September 2024 (05P)</p> <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - General information section rewritten. - The term 'mental retardation' is considered outdated and was updated to 'intellectual disability' where appropriate. - NM_ reference sequence for <i>ARID1B</i> gene has been updated according to MANE transcripts. - The ligation site and exon numbering for the ARID1B probes have been updated according to MANE Select transcript NM_001374828.1. - Added footnote "Δ More variable" for ARID1A 18899-L24494 and 18910-L32453 probes, and ARID1B 1893-L32281 probe to Table 1 and Table 2. - "Intron probe" footnote added to Table 2 for ARID1B 18939-L24534 probe. - List of References updated. - Various minor textual or layout changes. <p>Version A2-01 – 27 November 2020 (04P)</p> <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2). - Various minor textual or layout changes. - Ligation sites of the probes targeting the ARID1A and ARID1B genes updated according to new version of the NM_006015.6 and NM_001374820.1 reference sequences. - For uniformity, the chromosomal locations and bands in this document are now all based on hg18 (NCBI36).

More information: www.mrcholland.com ; www.mrcholland.eu	
	MRC Holland BV; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200

MRC Holland, SALSA, MLPA, digitalMLPA, Coffalyser.Net, Coffalyser digitalMLPA, and their logos are trademarks or registered trademarks of MRC Holland BV. All other brands and names herein are the property of their respective owners.