

Product Description

SALSA® MLPA® Probemix P075-B2 TCF4-FOXG1

To be used with the MLPA General Protocol.

Version B2

For complete product history see page 8.

Catalogue numbers:

- **P075-025R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 25 reactions.
- **P075-050R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 50 reactions.
- **P075-100R:** SALSA MLPA Probemix P075 TCF4-FOXG1, 100 reactions.

To be used in combination with a SALSA MLPA reagent kit and Coffalyser.Net data analysis software. MLPA reagent kits are either provided with FAM or Cy5.0 dye-labelled PCR primer, suitable for Applied Biosystems and Beckman/SCIEX capillary sequencers, respectively (see www.mrcholland.com).

Certificate of Analysis

Information regarding storage conditions, quality tests, and a sample electropherogram from the current sales lot is available at www.mrcholland.com.

Precautions and warnings

For professional use only. Always consult the most recent product description AND the MLPA General Protocol before use: www.mrcholland.com. It is the responsibility of the user to be aware of the latest scientific knowledge of the application before drawing any conclusions from findings generated with this product.

General information

The SALSA MLPA Probemix P075 TCF4-FOXG1 is a **research use only (RUO)** assay for the detection of deletions or duplications in the *TCF4* and *FOXG1* genes, which are associated with Pitt-Hopkins syndrome (PTHS).

PTHS is a rare disorder characterized by severe intellectual disability, pervasive developmental delay, 'atypical' autistic characteristics, and hyperventilation. PTHS is caused by heterozygous hypomorphic or null mutation or deletion of the transcription factor 4 (*TCF4*; *E2-2*; *ITF2*) gene on human chromosome 18. The *TCF4* gene is also a risk factor with highly significant linkage to schizophrenia, presumably via overexpression of the *TCF4* gene product in the central nervous system.

Another disorder that has some phenotypical overlap with PTHS is the congenital variant of Rett syndrome. This syndrome occurs almost exclusively in females and is characterized by hypotonia and intellectual disability from the very first months of life. It has been reported that the main cause of the congenital variant of Rett syndrome is mutations in the *FOXG1* gene, which encodes a transcription factor of the forkhead family on chromosome 14.

More information is available at <https://www.ncbi.nlm.nih.gov/books/NBK100240/>.

This SALSA MLPA probemix is not CE/FDA registered for use in diagnostic procedures. Purchase of this product includes a limited license for research purposes.

Gene structure and transcript variants:

Entrez Gene shows transcript variants of each gene: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=gene>
For NM_ mRNA reference sequences: <http://www.ncbi.nlm.nih.gov/sites/entrez?db=nucleotide>
Locus Reference Genomic (LRG) database: <http://www.lrg-sequence.org/>

Exon numbering

The *TCF4* and *FOXG1* exon numbering used in this P075-B2 TCF4-FOXG1 product description is the exon numbering from the NG_011716.2 and NG_009367.1 sequences, respectively. The exon numbering of the NM_ sequence that was used for determining a probe's ligation site does not always correspond to the exon numbering obtained from the LRG sequences. As changes to the databases can occur after release of this product description, the NM_ sequence and exon numbering may not be up-to-date.

Probemix content

The SALSA MLPA Probemix P075-B2 TCF4-FOXG1 contains 50 MLPA probes with amplification products between 130 and 485 nucleotides (nt). This includes 35 probes for the *TCF4* gene covering all exons, with additional probes for exons 1, 4, 5, 6, 9, 11, and 20, and four intronic probes and four probes located upstream of *TCF4*, most of them detecting exons of other transcript variants. Furthermore, five probes for the *FOXG1* are included, two probes upstream of the gene and three probes targeting exon 1. In addition, ten reference probes are included that detect autosomal chromosomal locations. Complete probe sequences and the identity of the genes detected by the reference probes are available online (www.mrcholland.com).

This probemix contains nine quality control fragments generating amplification products between 64 and 105 nt: four DNA Quantity fragments (Q-fragments), two DNA Denaturation fragments (D-fragments), one Benchmark fragment, and one chromosome X and one chromosome Y-specific fragment (see table below). More information on how to interpret observations on these control fragments can be found in the MLPA General Protocol and online at www.mrcholland.com.

Length (nt)	Name
64-70-76-82	Q-fragments (only visible with <100 ng sample DNA)
88-96	D-fragments (low signal indicates incomplete denaturation)
92	Benchmark fragment
100	X-fragment (X chromosome specific)
105	Y-fragment (Y chromosome specific)

MLPA technique

The principles of the MLPA technique (Schouten et al. 2002) are described in the MLPA General Protocol (www.mrcholland.com).

MLPA technique validation

Internal validation of the MLPA technique using 16 DNA samples from healthy individuals is required, in particular when using MLPA for the first time, or when changing the sample handling procedure, DNA extraction method or instruments used. This validation experiment should result in a standard deviation ≤ 0.10 for all probes over the experiment.

Required specimens

Extracted DNA free from impurities known to affect MLPA reactions. For more information please refer to the section on DNA sample treatment found in the MLPA General Protocol.

Reference samples

A sufficient number (≥ 3) of reference samples should be included in each MLPA experiment for data normalisation. All samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible. Reference samples should be derived from different unrelated individuals who are from families without a history of Pitt-Hopkins syndrome or Rett syndrome. More information regarding the selection and use of reference samples can be found in the MLPA General Protocol (www.mrcholland.com).

Positive control DNA samples

MRC Holland cannot provide positive DNA samples. Inclusion of a positive sample in each experiment is recommended. Coriell Institute (<https://catalog.coriell.org>) and Leibniz Institute DSMZ

(<https://www.dsmz.de/>) have diverse collections of biological resources which may be used as positive control DNA samples in your MLPA experiments. The quality of cell lines can change; therefore samples should be validated before use.

Data analysis

Coffalyser.Net software should be used for data analysis in combination with the appropriate lot-specific MLPA Coffalyser sheet. For both, the latest version should be used. Coffalyser.Net software is freely downloadable at www.mrcholland.com. Use of other non-proprietary software may lead to inconclusive or false results. For more details on MLPA quality control and data analysis, including normalisation, see the Coffalyser.Net Reference Manual.

Interpretation of results

The standard deviation of each individual probe over all the reference samples should be ≤ 0.10 and the final ratio (FR) of each individual reference probe in the patient samples should be between 0.80 and 1.20. When these criteria are fulfilled, the following cut-off values for the FR of the probes can be used to interpret MLPA

Copy number status	Final ratio (FR)
Normal	$0.80 < FR < 1.20$
Homozygous deletion	FR = 0
Heterozygous deletion	$0.40 < FR < 0.65$
Heterozygous duplication	$1.30 < FR < 1.65$
Heterozygous triplication/homozygous duplication	$1.75 < FR < 2.15$
Ambiguous copy number	All other values

Note: The term “dosage quotient”, used in older product description versions, has been replaced by “final ratio” to become consistent with the terminology of the Coffalyser.Net software. (Calculations, cut-offs and interpretation remain unchanged.) Please note that the Coffalyser.Net software also shows arbitrary borders as part of the statistical analysis of results obtained in an experiment. As such, arbitrary borders are different from the final ratio cut-off values shown here above.

- Arranging probes according to chromosomal location facilitates interpretation of the results and may reveal more subtle changes such as those observed in mosaic cases. Analysis of parental samples may be necessary for correct interpretation of complex results.
- False positive results: Please note that abnormalities detected by a single probe (or multiple consecutive probes) still have a considerable chance of being a false positive result. Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can be one cause. Incomplete DNA denaturation (e.g. due to salt contamination) can also lead to a decreased probe signal, in particular for probes located in or near a GC-rich region or in or near the *FOXG1* gene. The use of an additional purification step or an alternative DNA extraction method may resolve such cases. Additionally, contamination of DNA samples with cDNA or PCR amplicons of individual exons can lead to an increased probe signal (Varga et al. 2012). Analysis of an independently collected secondary DNA sample can exclude these kinds of contamination artefacts.
- Normal copy number variation in healthy individuals is described in the database of genomic variants: <http://dgv.tcag.ca/dgv/app/home>. Users should always consult the latest update of the database and scientific literature when interpreting their findings.
- Not all abnormalities detected by MLPA are pathogenic. In some genes, intragenic deletions are known that result in very mild or no disease (as described for *DMD* by Schwartz et al. 2007). For many genes, more than one transcript variant exists. Copy number changes of exons that are not present in all transcript variants may not have clinical significance. Duplications that include the first or last exon of a gene (e.g. exons 1-3) might not result in inactivation of that gene copy.
- Copy number changes detected by reference probes or flanking probes are unlikely to have any relation to the condition tested for.
- False results can be obtained if one or more peaks are off-scale. For example, a duplication of one or more exons can be obscured when peaks are off-scale, resulting in a false negative result. The risk on off-scale

peaks is higher when probemixes are used that contain a relatively low number of probes. Coffalyser.Net software warns for off-scale peaks while other software does not. If one or more peaks are off-scale, rerun the PCR products using either: a lower injection voltage or a shorter injection time, or a reduced amount of sample by diluting PCR products.

Limitations of the procedure

- In most populations, the major cause of genetic defects in the *TCF4* and *FOXP1* genes are small (point) mutations, most of which will not be detected by using SALSA MLPA Probemix P075 TCF4-FOXP1.
- MLPA cannot detect any changes that lie outside the target sequence of the probes and will not detect copy number neutral inversions or translocations. Even when MLPA did not detect any aberrations, the possibility remains that biological changes in that gene or chromosomal region *do* exist but remain undetected.
- Sequence changes (e.g. SNVs, point mutations) in the target sequence detected by a probe can cause false positive results. Mutations/SNVs (even when >20 nt from the probe ligation site) can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA.

Confirmation of results

Copy number changes detected by only a single probe always require confirmation by another method. An apparent deletion detected by a single probe can be due to e.g. a mutation/polymorphism that prevents ligation or destabilises the binding of probe oligonucleotides to the DNA sample. Sequence analysis can establish whether mutations or polymorphisms are present in the probe target sequence. The finding of a heterozygous mutation or polymorphism indicates that two different alleles of the sequence are present in the sample DNA and that a false positive MLPA result was obtained.

Copy number changes detected by more than one consecutive probe should be confirmed by another independent technique such as long range PCR, qPCR, array CGH or Southern blotting, whenever possible. Deletions/duplications of more than 50 kb in length can often be confirmed by FISH.

LOVD mutation database

<https://databases.lovd.nl/shared/genes/>. We strongly encourage users to deposit positive results in the Leiden Open Variation Database (LOVD). Recommendations for the nomenclature to describe deletions/duplications of one or more exons can be found on <http://varnomen.hgvs.org/>.

Please report copy number changes detected by the reference probes, false positive results due to SNVs and unusual results (e.g., a duplication of *TCF4* exons 13 and 15 but not exon 14) to MRC Holland: info@mrcholland.com.

Table 1. SALSA MLPA Probemix P075-B2 TCF4-FOXG1

Length (nt)	SALSA MLPA probe	Chromosomal position (hg18) ^a		
		Reference	TCF4	FOXG1
64-105	Control fragments – see table in probemix content section for more information			
130	Reference probe 00797-L13645	5q		
138	TCF4 probe 13335-L14761		Exon 4	
143	TCF4 probe 13334-L26942		Exon 14	
148	TCF4 probe 13337-L14763		Exon 18	
154	Reference probe 08375-L08229	15q		
160	TCF4 probe 13336-L14762		Exon 20	
166	TCF4 probe 17729-L21843		Upstream	
172	Reference probe 11007-L11678	4q		
178	TCF4 probe 16847-L19641		Exon 5	
184	TCF4 probe 17345-L20919		Exon 15	
190 «	FOXG1 probe 13756-L15243			Upstream
196	TCF4 probe 13333-L14759		Exon 6	
202	TCF4 probe 13339-L14765		Exon 8	
208 «	TCF4 probe 12506-L13556		Exon 1	
214	TCF4 probe 19805-L15889		Exon 4	
220	TCF4 probe 13342-L15879		Exon 16	
226 «	TCF4 probe 17730-L21844		Upstream	
232	Reference probe 09641-L26943	17q		
238	TCF4 probe 13326-L14752		Exon 10	
247 «	FOXG1 probe 17346-L20920			Exon 1
256	TCF4 probe 13324-L14750		Exon 6	
264	Reference probe 09265-L10877	10q		
274 «	FOXG1 probe 13755-L15242			Exon 1
283	TCF4 probe 13346-L14772		Exon 7	
292	TCF4 probe 13325-L14751		Exon 11	
301	TCF4 probe 17731-L21845		Upstream	
310 «	TCF4 probe 13341-L14767		Exon 1	
316 Ж Ø «	TCF4 probe 19601-SP0845-L26944		Intron 2	
322	TCF4 probe 13345-L26945		Exon 17	
328	Reference probe 08881-L08937	7q		
335 «	FOXG1 probe 16850-L27389			Upstream
342 Ø	TCF4 probe 19602-L27390		Intron 6	
350 «	TCF4 probe 16851-L27392		Exon 3	
358	TCF4 probe 13327-L27395		Exon 11	
364	Reference probe 12656-L19936	16q		
372	TCF4 probe 16852-L19646		Exon 19	
379	TCF4 probe 19603-L26947		Intron 5	
386	TCF4 probe 13329-L26946		Exon 12	
391	TCF4 probe 12522-L13572		Exon 20	
401	TCF4 probe 13343-L14769		Exon 5	
409	TCF4 probe 13344-L14770		Exon 9	
418	Reference probe 06876-L05967	3p		
427 «	TCF4 probe 16853-L19647		Exon 2	
436	TCF4 probe 13340-L14766		Exon 9	
445	Reference probe 10667-L11249	6p		
454	TCF4 probe 19604-L26951		Intron 8	
461	TCF4 probe 13348-L26950		Exon 13	
468 «	FOXG1 probe 13754-L26949			Exon 1
476 Ж	TCF4 probe 17732-SP0544-L26948		Upstream	
485	Reference probe 13594-L22376	19p		

^a See section Exon numbering on page 2 for more information.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Table 2. P075-B2 probes arranged according to chromosomal location

Table 2a. *TCF4*

Length (nt)	SALSA MLPA probe	<i>TCF4</i> exon ^a	Ligation site NM_001083962.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		<i>start codon</i>	176-178 (<i>Exon 2</i>)		
476 Ж	17732-SP0544-L26948	Upstream	NM_001243226.3; 97-98 and 130-131	TAGCACAGGTGC-33 nt spanning oligo-CAAATTACTCAA	0.1 kb
301	17731-L21845	Upstream	NM_001243226.3; 235-236	GCTACTCGAGCT-TCTCCAGGAGGC	4.0 kb
166	17729-L21843	Upstream	NM_001243226.3; 315 nt before exon 2	CTAAGGCACTCA-CCTACTCTCAGA	41.9 kb
226 «	17730-L21844	Upstream	NM_001243227.2; 46-45, reverse	TCCCTGAAAGAT-ACATTGTAATCC	1.2 kb
208 «	12506-L13556	Exon 1	338 nt before exon 1	CCGAGGGATGCA-ACGGGCAAAAAC	0.3 kb
310 «	13341-L14767	Exon 1	35 nt before exon 1, reverse	CTCTTAACACCA-ACTCTCTTCTCC	1.2 kb
427 «	16853-L19647	Exon 2	226-225, reverse	TCCAGTAAATCA-CTCAGCTCTTTG	1.1 kb
316 Ж ∅ «	19601-SP0845-L26944	Intron 2	698 nt and 668 nt before exon 3	CCCTAGGCAGGC-30 nt spanning oligo-CTTTCTCCATTC	0.7 kb
350 «	16851-L27392	Exon 3	305-306	TGGCAAGTGGAC-ATTTTACTGGCT	121.0 kb
214	19805-L15889	Exon 4	147 nt before exon 4, reverse	AGTGGCTTCTGA-CCCATCTACTTA	0.2 kb
138	13335-L14761	Exon 4	14 nt after exon 4, reverse	TGGGAGAAAAGA-TTAGATATACTT	2.8 kb
401	13343-L14769	Exon 5	130 nt before exon 5	GATTCCTTCTAG-TGAAGTTCCAGG	0.2 kb
178	16847-L19641	Exon 5	446-447	CACATGACAATC-TCTCTCCACCTT	38.8 kb
379	19603-L26947	Intron 5	18.8 kb before exon 6 (NM_001243233.2; 20-21)	GCCACAACAGTT-TATTCATCCACA	18.7 kb
196	13333-L14759	Exon 6	100 nt before exon 6	CCTACTTTACGT-ATGTAAACATCG	0.1 kb
256	13324-L14750	Exon 6	512-513	ACTCATCTTATG-GGAGAGAATCAA	1.3 kb
342 ∅	19602-L27390	Intron 6	1.3 kb after exon 6	ACAGTGCTTGGT-TAAGAGCTCCTG	51.2 kb
283	13346-L14772	Exon 7	636-635, reverse	TTCGGGGATTAT-TGCTAGAATACT	0.5 kb
202	13339-L14765	Exon 8	703-702, reverse	AAACCTGGAGGA-ACTTTTTCGAACT	28.7 kb
454	19604-L26951	Intron 8	28.7 kb after exon 8 (NM_001243234.2; 119-120)	ACTGCGCATACA-CAATCCCGGGCA	42.1 kb
436	13340-L14766	Exon 9	740-741	ATGCTCCATCAG-CAAGCACTGCCG	0.1 kb
409	13344-L14770	Exon 9	8 nt after exon 9	CAAGGTAAGATG-CTGCTGCTTCTG	3.9 kb
238	13326-L14752	Exon 10	919-920	TCTTCTCATATT-CCACAGTCCAGC	5.8 kb
292	13325-L14751	Exon 11	1020-1021	TCCGATGTCCAC-TTTCATCGTAG	0.1 kb
358	13327-L27395	Exon 11	18 nt after exon 11	CACAGAAATGCC-AATTCTGATACC	8.3 kb
386	13329-L26946	Exon 12	1162-1161, reverse	TTCTCACCGAA-GCAAGTGCTTTC	1.6 kb
461	13348-L26950	Exon 13	17 nt after exon 13	GTATTTCAAATC-CCATTTTCATCAT	2.5 kb
143	13334-L26942	Exon 14	1247-1248	TAATATCAGCAG-GCACAGCTGTTT	2.8 kb
184	17345-L20919	Exon 15	1420-1419, reverse	TGCATGTCCCA-TGACCACCAGGC	20.0 kb

Length (nt)	SALSA MLPA probe	TCF4 exon ^a	Ligation site NM_001083962.2	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
220	13342-L15879	Exon 16	1609-1610	CCACAGCTTCCT-GTCCAGTCTGCG	1.9 kb
322	13345-L26945	Exon 17	60 nt before exon 17	GCAGCCTTGCAA-TCTGGTGTGCAG	3.7 kb
148	13337-L14763	Exon 18	1827-1828	GTCCCACAGCAA-TAATGACGATGA	0.8 kb
372	16852-L19646	Exon 19	2156-2157	CACACCCTGGAA-TGGGAGACGCAT	0.5 kb
160	13336-L14762	Exon 20	2467-2468	ACAGGCTGAGAC-ACAGCCCAGAGA	0.6 kb
391	12522-L13572	Exon 20	3028-3029	CCTGTAGTGCCA-ACTCTGCTTCCA	
		stop codon	2189-2191 (Exon 19)		

Table 2b. *FOXG1*

Length (nt)	SALSA MLPA probe	FOXG1 exon ^a	Ligation site NM_005249.5	Partial sequence ^b (24 nt adjacent to ligation site)	Distance to next probe
		start codon	494-496 (Exon 1)		
335 «	16850-L27389	Upstream	559 nt before exon 1	AAATGCCAGACA-CTGGCCTGCAAG	0.2 kb
190 «	13756-L15243	Upstream	349 nt before exon 1	GAGGAAGCCGGA-AATGTGAGCTAT	1.9 kb
247 «	17346-L20920	Exon 1	1508-1509	CCAGCCACCCCA-TGCCCTACAGCT	0.4 kb
274 «	13755-L15242	Exon 1	1907-1906, reverse	GAAATAATCAGA-CAGTCCCCCAGA	0.2 kb
468 «	13754-L26949	Exon 1	2101-2102	TCTAGGGTTGTT-TATTATTCTAAC	
		stop codon	1961-1963 (Exon 1)		

^a See section Exon numbering on page 2 for more information.

^b Only partial probe sequences are shown. Complete probe sequences are available at www.mrcholland.com. Please notify us of any mistakes: info@mrcholland.com.

« Probe located in or near a GC-rich region. A low signal can be caused by salt contamination in the DNA sample leading to incomplete DNA denaturation, especially of GC-rich regions.

Ж This probe consists of three parts and has two ligation sites. A low signal of this probe can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples.

∅ Intron probe. Only included to help determine the extent of a deletion/duplication. Copy number alterations of only this probe are of unknown clinical significance.

SNVs located in the target sequence of a probe can influence probe hybridization and/or probe ligation. Single probe aberration(s) must be confirmed by another method.

Complete probe sequences are available at www.mrcholland.com.

Related SALSA MLPA probemixes

ME028 Prader Willi/Angelman	Contains probes for different genes and regions related to Prader-Willi and Angelman syndrome.
P015 MECP2	Contains probes for the <i>MECP2</i> gene and Xq28 region related to RETT syndrome.
P169 Hirschsprung-1	Contains probes for the <i>ZEB2</i> gene (Mowat-Wilson) and other genes related to Hirschsprung disease.
P189 CDKL5	Contains probes for different genes and regions related to RETT-like syndrome.
P336 UBE3A	Contains probes for the <i>UBE3A</i> , <i>MTHFR</i> , and <i>GABRB3</i> genes which are related to Angelman syndrome.
P379 NRXN1	Contains probes for the <i>NRXN1</i> gene related to Pitt-Hopkins-like syndrome 2.

References

- Schouten JP et al. (2002). Relative quantification of 40 nucleic acid sequences by multiplex ligation-dependent probe amplification. *Nucleic Acids Res.* 30:e57.

- Schwartz M et al. (2007). Deletion of exon 16 of the dystrophin gene is not associated with disease. *Hum Mutat.* 28:205.
- Varga RE et al. (2012). MLPA-based evidence for sequence gain: pitfalls in confirmation and necessity for exclusion of false positives. *Anal Biochem.* 421:799-801.

Selected publications using SALSA MLPA Probemix P075 TCF4-FOXG1

- Bedeschi MF et al. (2017). Impairment of different protein domains causes variable clinical presentation within Pitt-Hopkins syndrome and suggests intragenic molecular syndromology of TCF4. *Eur J Med Genet* 60.11: 565-571.
- Cellini E et al. (2016). The hyperkinetic movement disorder of FOXG1-related epileptic–dyskinetic encephalopathy. *Dev Med Child Neurol* 58.1: 93-97.
- Currò A et al. (2021). CDKL5 mutations may mimic Pitt-Hopkins syndrome phenotype. *Eur J Med Genet*, 64(1), 104102.
- De Bruyn C et al. (2014). Thin genu of the corpus callosum points to mutation in FOXG1 in a child with acquired microcephaly, trigonocephaly, and intellectual developmental disorder: a case report and review of literature. *Eur J Paediatr Neurol*, 18(3), 420-426.
- Hellwig M et al. (2019). TCF4 (E2-2) harbors tumor suppressive functions in SHH medulloblastoma. *Acta neuropathol*, 137(4), 657-673.
- Kumakura A et al. (2014). A haploinsufficiency of FOXG1 identified in a boy with congenital variant of Rett syndrome. *Brain Dev*, 36(8), 725-729.
- Tripon F et al. (2020). Pitt-Hopkins Syndrome: Clinical and Molecular Findings of a 5-Year-Old Patient. *Genes*, 11(6), 596.
- Vidal S et al. (2017). The utility of Next Generation Sequencing for molecular diagnostics in Rett syndrome. *Sci Rep*, 7(1), 1-11.

P075 product history	
Version	Modification
B2	Two reference probes have been replaced and one probe length has been adjusted.
B1	One <i>FOXG1</i> probe is removed.
A1	First release.

Implemented changes in the product description
<p>Version B2-03 – 25 November 2024 (04P)</p> <ul style="list-style-type: none"> - The term ‘mental retardation’ is considered outdated and was updated to ‘intellectual disability’ where appropriate.
<p>Version B2-02 – 10 November 2021 (04P)</p> <ul style="list-style-type: none"> - Product description rewritten and adapted to a new template. - Ligation sites of the probes targeting the <i>TCF4</i> gene updated according to new version of the NM_ sequences. - Small changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.
<p>Version B2-01 – 21 February 2018 (01P)</p> <ul style="list-style-type: none"> - Product description restructured and adapted to a new template. - Product description adapted to a new product version (version number changed, changes in Table 1 and Table 2, new picture included).
<p>Version 08 (53)</p> <ul style="list-style-type: none"> - For the <i>TCF4</i> and <i>FOXG1</i> gene exon numbers and ligation sites in table 1 and 2 have been adjusted according to NCBI Map Viewer. - Product description adapted to a new product version (version number changed, lot number added, new picture included). - Various textual changes on page 1 and 2. - Changes of probe lengths in Table 1 and 2 in order to better reflect the true lengths of the amplification products.

Version 07 (48)

- Electropherogram pictures using the new MLPA buffer (introduced in December 2012) added.


Version 06 (48)

- Exon numbering of the *TCF4* gene has been changed due to the identification of new exons
- Various textual changes on page 1 related to the exon numbering change.

Version 05 (48)

- Remark on RefSeqGene standard and transcript variant added below Table 2.
- Various minor textual changes on page 1.
- Ligation sites of the probes targeting the *TCF4* and *FOXP1* genes updated according to new version of the NM_reference sequence.
- Small correction of chromosomal locations in Table 1 and 2.

More information: www.mrcholland.com; www.mrcholland.eu

	MRC Holland bv; Willem Schoutenstraat 1 1057 DL, Amsterdam, The Netherlands
E-mail	info@mrcholland.com (information & technical questions) order@mrcholland.com (orders)
Phone	+31 888 657 200