

Instructions for Use

SALSA® MLPA® Probemix P051 Parkinson mix 1 and SALSA® MLPA® Probemix P052 Parkinson mix 2



See also the MLPA General Protocol, and the product descriptions of the SALSA® MLPA® Reagent Kit, SALSA® Binning DNA SD067, and the Coffalyser.Net Reference Manual.

Visit the SALSA® MLPA® Probemix P051 Parkinson mix 1 and SALSA® MLPA® Probemix P052 Parkinson mix 2 product pages on our website to find Certificates of Analysis and a list of related products.

| | |
|--------------------------|---|
| Product Name | SALSA® MLPA® Probemix P051 Parkinson mix 1 |
| Version | D2 |
| Catalogue numbers | P051-025R (25 reactions) P051-050R (50 reactions) P051-100R (100 reactions) |
| Basic UDI-DI: | 872021148P0515M |
| Ingredients | Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA |


| | |
|--------------------------|---|
| Product Name | SALSA® MLPA® Probemix P052 Parkinson mix 2 |
| Version | D2 |
| Catalogue numbers | P052-025R (25 reactions) P052-050R (50 reactions) P052-100R (100 reactions) |
| Basic UDI-DI: | 872021148P0525P |
| Ingredients | Synthetic oligonucleotides, oligonucleotides purified from bacteria, Tris-HCl, EDTA |

| Additional Test Components | Catalogue numbers |
|--|--|
| SALSA® MLPA® Reagent Kit | EK1-FAM EK1-CY5 EK5-FAM EK5-CY5 EK20-FAM |
| SALSA® Binning DNA SD067 | SD067 |


Storage and Shelf Life

| | | |
|------------------------|---|---|
| Recommended conditions |  |  |
|------------------------|---|---|

A shelf life of until the expiry date is guaranteed, also after opening when stored in the original packaging under recommended conditions. For the exact expiry date, see the label on the vial. These products should not be exposed to more than 25 freeze-thaw cycles. Do not use the products if the packaging is damaged or opened. Leave chemicals in original containers. Waste material must be disposed of in accordance with the national and local regulations.

| Regulatory Status | |
|-------------------|---|
| IVD | EUROPE  2797 ISRAEL |
| RUO | ALL OTHER COUNTRIES |

| Label Symbols | | | |
|---------------|---------------------|------------|-------------------|
| IVD | In Vitro Diagnostic | RUO | Research Use Only |

| More Information: | |
|---|--|
| www.mrcholland.com | |
|  | MRC Holland BV; Willem Schoutenstraat 1 1057 DL, Amsterdam, the Netherlands |
| E-mail | info@mrcholland.com (information & technical questions); order@mrcholland.com (orders) |
| Phone | +31 888 657 200 |

Any serious incident that has occurred in relation to these products should be reported to MRC Holland and the competent authority of the Member State in which the user and/or the patient is located.

Changes in these Product Versions

P051 version D2. As compared to version D1, two probes have a change in length and sequence.

P052 version D2. As compared to version D1, one probe has a change in length and sequence and one probe has a change in sequence.

1. Intended Purpose

The SALSA MLPA Probemix P051 Parkinson mix 1 and SALSA MLPA Probemix P052 Parkinson mix 2 are in vitro diagnostic (IVD)¹ or research use only (RUO) semi-quantitative manual assays² for the detection of duplications and triplications in the *SNCA* gene (P051), deletions and duplications in the *PARK2* gene, and the presence of one point mutation, G2019S in the *LRRK2* gene (P051 and P052), in genomic DNA isolated from human peripheral whole blood specimens. P051 Parkinson mix 1 and P052 Parkinson mix 2 are intended to confirm a potential cause for early-onset (*PARK2* deletions/duplication, *SNCA* triplications, *LRRK2* G2019S mutation) and late-onset (*SNCA* duplications, *LRRK2* G2019S mutations) Parkinson’s disease and for molecular genetic testing of at-risk family members³.

Copy number variations (CNVs) and the *LRRK2* point mutation, G2019S, detected with P051 Parkinson mix 1 and P052 Parkinson mix 2 should be confirmed with a different technique. In particular, CNVs detected by only a single probe always require confirmation by another method. Most defects in *SNCA* and *PARK2* are point mutations, which will not be detected by MLPA, with the exception of the aforementioned *LRRK2* G2019S mutation. It is therefore recommended to use this assay in combination with sequence analysis.

Assay results are intended to be used in conjunction with other clinical and diagnostic findings, consistent with professional standards of practice, including confirmation by alternative methods, clinical genetic evaluation, and counselling, as appropriate. The results of this test should be interpreted by a clinical molecular geneticist or equivalent.

These devices are not intended to be used for standalone diagnostic purposes, pre-implantation or prenatal testing, population screening, or for the detection of, or screening for, acquired or somatic genetic aberrations.

¹ Please note that these probemixes are for in vitro diagnostic (IVD) use in the countries specified on page 1 of this product description. In all other countries, the products are research use only (RUO).

² To be used in combination with a SALSA MLPA Reagent Kit, Coffalyser.Net analysis software, and SALSA Binning DNA SD067.

³ Certain probes included in P051 and P052 targeting other genes may only be used in a research setting. The following table summarises which probes are for IVD or exclusively restricted to RUO use:

| | IVD targets | RUO targets only |
|-------------|---|---|
| P051 | CNVs: <i>SNCA</i> ; <i>PARK2</i> Mutation: <i>LRRK2</i> G2019S | CNVs: <i>PARK7</i> , <i>ATP13A2</i> , <i>PINK</i> Mutation: <i>SNCA</i> A30P |
| P052 | CNVs: <i>PARK2</i> Mutation: <i>LRRK2</i> G2019S | CNVs: <i>UCHL1</i> , <i>LRRK2</i> , <i>GCH1</i> , <i>ATP13A2</i> , <i>CAV1/2</i> |

2. Sample Requirements

| | |
|-------------------|--|
| Specimen | 50-250 ng purified human genomic DNA, free from heparin, dissolved in 5 µl TE _{0.1} buffer, pH 8.0-8.5 |
| Collection Method | Standard methods |
| Extraction Method | Methods tested by MRC Holland: <ul style="list-style-type: none"> • QIAGEN Autopure LS (automated) and QIAamp DNA mini/midi/maxi kit (manual) • Promega Wizard Genomic DNA Purification Kit (manual) • Salting out (manual) |

| Sample Types | | |
|---------------------------------------|--|------------------------------|
| Test Sample | <ul style="list-style-type: none"> • Provided by user | |
| Reference Samples (Required) | <ul style="list-style-type: none"> • Provided by user • Extraction method, tissue type, DNA concentration (and) treatment as similar as possible in all test and reference samples. • Have a normal copy number and ≤ 0.10 standard deviation for all probes except for mutation-specific probes. • At least three* independent reference samples required in each experiment for proper data normalisation. Derived from unrelated individuals from families without a history of Parkinson’s Disease (PD). | |
| No-DNA Control (Preferably) | <ul style="list-style-type: none"> • Provided by user • TE_{0.1} buffer instead of DNA • To check for DNA contamination | |
| Binning Sample (Initial Experiment) | <ul style="list-style-type: none"> • SALSA® Binning DNA SD067, provided by MRC Holland • Required in initial experiment to determine suitable bin set • Should never be used as a reference sample | |
| Positive Control Samples (Preferably) | <ul style="list-style-type: none"> • Provided by user, or | |
| | <table border="1"> <tr> <td>Available from third parties</td> <td>See the table of positive samples on the probemixes product pages on our website.</td> </tr> </table> | Available from third parties |
| Available from third parties | See the table of positive samples on the probemixes product pages on our website. | |
| Validation Samples (Required) | <ul style="list-style-type: none"> • In the validation experiments of this probemix, the peaks of the mutation-specific probes are expected to be absent in the majority of samples from healthy individuals. | |

*When testing >21 samples, include one extra reference for each 7 test samples.

3. Test Procedure

See the [MLPA General Protocol](#).

4. Quality Control, Data Analysis, and Troubleshooting

| Quality Control Fragments in the Probemixes | |
|---|--|
| Length (nt) | Function |
| 64-70-76-82 | DNA Quantity control fragments |
| 88-96 | DNA Denaturation control fragments |
| 92 | Benchmark fragment |
| 100 | Chromosome X presence control fragment |
| 105 | Chromosome Y presence control fragment |

[Coffalyser.Net](#) should be used for data analysis in combination with the appropriate product and lot-specific Coffalyser sheet. See the [Coffalyser.Net Reference Manual](#) for details on data analysis and quality control.

For troubleshooting help, see the additional resources offered on our [support portal](#).

5. Interpretation of Results

Determining Typical Values in Normal and Affected Populations

The typical final ratio (FR) values stated in the copy number tables were determined in a validation study with samples containing abnormal copy numbers. The standard deviation of each individual probe over all the reference samples was ≤ 0.10 .

6. Performance Characteristics

| Study | Description |
|---|---|
| Expected values for copy numbers in normal and affected populations | To determine the expected values in normal and affected populations a study was conducted on over 1500 MLPA reactions using samples with and without abnormal copy numbers. When the standard deviation of each individual probe over all the reference samples is ≤ 0.10 , the ranges stated in the copy number table can be used. Cut-off values for copy number determination were verified with P051 in 46 samples from healthy individuals with normal copy number and eight samples with known CNVs and with P052 in 43 samples from healthy individuals with normal copy number and five samples with known CNVs. The expected FRs for the corresponding copy number were found in all samples tested. |
| Expected values for point mutation detection in normal and affected populations | The mutation-specific probe will only generate a signal when the <i>LRRK2</i> sG2019S (196 nt for P051 and 172 nt for P052) mutation is present. Please note that background signals of the mutation-specific probes can be expected above the threshold in some cases. Users should always compare the relative peak height of the mutation-specific probe in mutation-positive samples to the relative peak height in reference samples. A clear signal (at least 10% of the median peak height of all reference probes in that sample) indicates that the mutation is present. It is not possible to determine the copy number of mutation-specific probes. The expected value for mutation-specific probes was verified with P051 and P052 using one positive sample for the <i>LRRK2</i> G2019S mutation and 46 and 43 samples from healthy individuals without the mutation, respectively, and the expected mutation status was found in each case. |
| Limit of detection | A study that evaluated the acceptable minimum and maximum amount of sample DNA revealed that the use of 50-250 ng of human DNA is the recommended input. The use of insufficient or too much sample DNA can affect performance. These lower and higher limits of detection were verified using P051 on one sample with known CNVs and using P052 on four samples with known CNVs/mutations and on one sample without any mutation and expected results were obtained using both the lower and upper input amount of DNA. |
| Interfering substances | Impurities in the DNA sample can affect the MLPA reaction. To minimise this effect, see Sample quality section under Precautions and warnings of the MLPA General Protocol. A study using P051 and P052 was performed to assess the potential for interference of endogenous and exogenous substances on genomic DNA from samples with known CNVs/mutations. For most probes, |

Expected Results of Reference Probes

| Final Ratio (FR) | Copy Number | Description |
|------------------|-------------|-------------|
| 0.80-1.20 | 2 | Normal |

Typical Results of Probes Targeting Two Copies (e.g. *SNCA*, *PARK2*)

| Final Ratio (FR) | Copy Number | Description |
|--------------------|-------------|---|
| 0 | 0 | Homozygous deletion |
| 0.40 – 0.65 | 1 | Heterozygous deletion |
| 0.80 – 1.20 | 2 | Normal |
| 1.30 – 1.65 | 3 | Heterozygous duplication |
| 1.75 – 2.15 | 4 | Homozygous duplication or Heterozygous triplication |
| All other values | - | Ambiguous |

The tables illustrate the relationship between final probe ratio and corresponding copy number. Test results are expected to centre around these values. Ambiguous values can indicate a technical problem, but may also reflect a biological cause such as mosaicism or a SNV influencing a single probe. It is important to use Coffalyser.Net to determine the significance of values found.

Possible Results of Mutation-Specific Probes

| Signal Strength | Mutation Status |
|---|---|
| $\geq 10\%$ median peak height reference probes | Mutation is detected (expected only in positive samples) |
| $< 10\%$ median peak height reference probes | Mutation is not detected (expected in most samples from healthy individuals) |

| | <p>expected FRs (FRs within the expected cut-off category) were obtained even in the presence of potential interferents at concentrations shown in the table below.</p> <table border="1" data-bbox="443 219 1437 884"> <thead> <tr> <th>Interferent</th> <th>Source</th> <th>Testing Concentration</th> <th>P051</th> <th>Results* P052</th> </tr> </thead> <tbody> <tr> <td>EDTA</td> <td>Exogenous – specimen collection tubes</td> <td>1.5 mM</td> <td>Copy number: Ambiguous FR for 4/108 probes Mutation: Expected status for 6/6 probes</td> <td>Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe</td> </tr> <tr> <td>NaCl</td> <td>Exogenous – DNA extraction</td> <td>40 mM</td> <td>Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes</td> <td>Copy number: Ambiguous FR for 2/50 probes Mutation: Expected status for 1/1 probe</td> </tr> <tr> <td>Fe³⁺ (FeCl₃)</td> <td>Exogenous – DNA extraction</td> <td>1 µM</td> <td>Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes</td> <td>Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe</td> </tr> <tr> <td>Heparin</td> <td>Exogenous – specimen collection tubes</td> <td>0.02 U/mL</td> <td>Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes</td> <td>Copy number: Ambiguous FR for 1/50 probes Mutation: Expected status for 1/1 probe</td> </tr> <tr> <td>Haemoglobin</td> <td>Endogenous – blood sample</td> <td>0.02 µg/µl</td> <td>Copy number: Expected FR for 17/108 probes Mutation: Expected status for 6/6 probes</td> <td>Copy number: Expected FR for 25/50 probes Mutation: Expected status for 1/1 probe</td> </tr> </tbody> </table> <p>* Results are summarised for all probes across all six and five samples tested for P051 and P052, respectively.</p> <p>Exogenous interfering substances (EDTA, heparin, salts (NaCl), and FeCl₃) were tested and shown to have a mild effect, leading to, at the most, ambiguous ratios and potential delayed results. Haemoglobin had the largest effect on the FRs, in particular for copy number determination.</p> <p>Additionally, Coffalyser.Net issues warnings for the samples in which the interferents showed an effect, as well as lowered quality scores.</p> <p>To minimise variability across samples, all samples tested, including reference DNA samples, should be derived from the same tissue type, handled using the same procedure, and prepared using the same DNA extraction method when possible.</p> | Interferent | Source | Testing Concentration | P051 | Results* P052 | EDTA | Exogenous – specimen collection tubes | 1.5 mM | Copy number: Ambiguous FR for 4/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe | NaCl | Exogenous – DNA extraction | 40 mM | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Ambiguous FR for 2/50 probes Mutation: Expected status for 1/1 probe | Fe ³⁺ (FeCl ₃) | Exogenous – DNA extraction | 1 µM | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe | Heparin | Exogenous – specimen collection tubes | 0.02 U/mL | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Ambiguous FR for 1/50 probes Mutation: Expected status for 1/1 probe | Haemoglobin | Endogenous – blood sample | 0.02 µg/µl | Copy number: Expected FR for 17/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 25/50 probes Mutation: Expected status for 1/1 probe |
|---------------------------------------|--|-----------------------|--|--|------|------------------|------|---------------------------------------|--------|--|--|------|----------------------------|-------|--|--|---------------------------------------|----------------------------|------|--|--|---------|---------------------------------------|-----------|--|--|-------------|---------------------------|------------|--|--|
| Interferent | Source | Testing Concentration | P051 | Results* P052 | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| EDTA | Exogenous – specimen collection tubes | 1.5 mM | Copy number: Ambiguous FR for 4/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| NaCl | Exogenous – DNA extraction | 40 mM | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Ambiguous FR for 2/50 probes Mutation: Expected status for 1/1 probe | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Fe ³⁺ (FeCl ₃) | Exogenous – DNA extraction | 1 µM | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 50/50 probes Mutation: Expected status for 1/1 probe | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Heparin | Exogenous – specimen collection tubes | 0.02 U/mL | Copy number: Ambiguous FR for 7/108 probes Mutation: Expected status for 6/6 probes | Copy number: Ambiguous FR for 1/50 probes Mutation: Expected status for 1/1 probe | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Haemoglobin | Endogenous – blood sample | 0.02 µg/µl | Copy number: Expected FR for 17/108 probes Mutation: Expected status for 6/6 probes | Copy number: Expected FR for 25/50 probes Mutation: Expected status for 1/1 probe | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Cross-reactivity | <p>Cross-reactivity is the potential for probes to bind to homologous regions (e.g. pseudogenes) or other cross-reactive sequences. Quality tests were carried out to determine whether probes are specific to their target sequence and all probes met the quality criteria for specificity.</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Accuracy | <p>Results of accuracy are derived from trueness and precision studies. Trueness: previously genotyped samples were tested and found to have the expected results. Precision: results are not affected by operator, day, or laboratory site. For trueness, ten previously genotyped samples were tested using P051, and six using P052, and found to have the expected results. Assay precision was tested by repeatedly testing samples with known copy number/mutations over multiple days, and by multiple operators. Overall, 99% correct calls (1610/1620 probes) for CNVs and 100% (90/90 probes) correct calls for mutation status were obtained throughout the precision experiments for P051. For P052, 99% correct calls (820/825 probes) for CNVs and 100% correct calls (75/75 probes) for mutation status were obtained throughout the precision experiments.</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| Clinical validity* | <ul style="list-style-type: none"> • <i>SNCA</i> (P051 only): less than 1% of Parkinson’s disease is caused by duplications and triplications in <i>SNCA</i> [1]. • <i>PARK2</i>: 2% of early onset and idiopathic, late-onset Parkinson’s disease is caused by deletions and duplications in <i>PARK2</i> [2]. • <i>LRRK2</i>: the frequency of <i>LRRK2</i> G2019S mutation varies across populations [3-5]. <p>*(Based on a 2000-2024 literature review)</p> | | | | | | | | | | | | | | | | | | | | | | | | | | | | | | |

Summary of Safety and Performance (SSP)

The SSP is available in the European database on medical devices (Eudamed), <https://ec.europa.eu/tools/eudamed>, or upon request.

Content – Probe Details Sorted by Chromosomal Position

SALSA MLPA Probemix P051-D2 Parkinson mix 1

| Chr. position | Target | Exon | Mutation | Distance to next probe | Length (nt) | Probe number | Warnings |
|---------------|----------------|-------------------|----------|------------------------|-------------|---------------------|----------|
| 1p36.23 | <i>TNFRSF9</i> | | | 20.6 kb | 136 | 20271-L27994 | - f |
| 1p36.23 | <i>PARK7</i> | Upstream (Exon 1) | | 0.3 kb | 429 | 03690-L27587 | ∅ |
| 1p36.23 | <i>PARK7</i> | Exon 1 | | 1.1 kb | 350 | 20279-L27649 | |
| 1p36.23 | <i>PARK7</i> | Exon 2 | | 2.6 kb | 272 | 20274-L27644 | |
| 1p36.23 | <i>PARK7</i> | Exon 3 | | 4.0 kb | 413 | 20101-L27586 | |
| 1p36.23 | <i>PARK7</i> | Exon 4 | | 1.6 kb | 319 | 20277-L27647 | |
| 1p36.23 | <i>PARK7</i> | Exon 5 | | 6.8 kb | 370 | 20254-L27588 | |
| 1p36.23 | <i>PARK7</i> | Exon 6 | | 7.4 kb | 245 | 20273-L27643 | |
| 1p36.23 | <i>PARK7</i> | Exon 7 | | 9.2 Mb | 457 | 02189-L27590 | |
| 1p36.13 | <i>ATP13A2</i> | Exon 9 | | 5.4 kb | 178 | 11715-L27546 | |
| 1p36.13 | <i>ATP13A2</i> | Exon 2 | | 3.6 Mb | 215 | 11716-L28013 | |
| 1p36.12 | <i>PINK1</i> | Exon 1 | | 4.1 kb | 149 | 21469-L30271 | Δ |
| 1p36.12 | <i>PINK1</i> | Exon 2 | | 2.1 kb | 172 | 03692-L27545 | |
| 1p36.12 | <i>PINK1</i> | Exon 3 | | 4.5 kb | 209 | 12067-L28012 | |
| 1p36.12 | <i>PINK1</i> | Exon 4 | | 1.2 kb | 405 | 20280-L27650 | |
| 1p36.12 | <i>PINK1</i> | Exon 5 | | 2.9 kb | 143 | 20270-L27992 | |
| 1p36.12 | <i>PINK1</i> | Exon 6 | | 0.5 kb | 310 | 20276-L27646 | |
| 1p36.12 | <i>PINK1</i> | Exon 7 | | 2.0 kb | 229 | 03698-L03154 | |
| 1p36.12 | <i>PINK1</i> | Exon 8 | | | 469 | 03697-L27591 | |
| 4q22.1 | <i>SNCA</i> | Exon 6 | | 2.6 kb | 184 | 02169-L28011 | |
| 4q22.1 | <i>SNCA</i> | Exon 5 | | 93.1 kb | 486 | 03689-L27592 | |
| 4q22.1 | <i>SNCA</i> | Exon 4 | | 5.8 kb | 166 | 02168-L27544 | |
| 4q22.1 | <i>SNCA</i> | Exon 3 | | 7.4 kb | 253 | 04616-L27552 | |
| 4q22.1 | <i>SNCA</i> | Exon 2 | A30P | 0.1 kb | 154 | 02166-L27543 | § |
| 4q22.1 | <i>SNCA</i> | Exon 2 | | 1.3 kb | 279 | 20255-L03103 | |
| 4q22.1 | <i>SNCA</i> | Exon 1 | | | 450 | 04096-L27589 | |
| 6q25.3 | <i>LPA</i> | | | 817.2 kb | 265 | 20224-L27548 | - † |
| 6q26 | <i>PARK2</i> | Exon 12 | | 10.4 kb | 395 | 02184-L27585 | |
| 6q26 | <i>PARK2</i> | Exon 11 | | 26.7 kb | 377 | 02183-L27896 | |
| 6q26 | <i>PARK2</i> | Exon 10 | | 162.0 kb | 359 | 02182-L27556 | |
| 6q26 | <i>PARK2</i> | Exon 9 | | 20.5 kb | 325 | 02181-L27555 | |
| 6q26 | <i>PARK2</i> | Exon 8 | | 216.5 kb | 302 | 02180-L27553 | |
| 6q26 | <i>PARK2</i> | Exon 7 | | 187.5 kb | 494 | 20283-L27895 | |
| 6q26 | <i>PARK2</i> | Exon 6 | | 80.8 kb | 343 | 20278-L27648 | |
| 6q26 | <i>PARK2</i> | Exon 5 | | 147.0 kb | 423 | 20281-L27651 | |
| 6q26 | <i>PARK2</i> | Exon 4 | | 61.4 kb | 202 | 20272-SP0951-L27900 | Ж |
| 6q26 | <i>PARK2</i> | Exon 3 | | 180.9 kb | 477 | 20282-L27652 | |
| 6q26 | <i>PARK2</i> | Exon 2 | | 284.3 kb | 287 | 02174-L27554 | |
| 6q26 | <i>PARK2</i> | Exon 1 | | | 237 | 20225-L24881 | |
| 12q12 | <i>LRRK2</i> | Exon 41 | G2019S | | 196 | 04575-L27549 | § |
| 2p | Reference | | | | 500 | 19555-L27674 | |
| 2q | Reference | | | | 335 | 18737-L27897 | |
| 3p | Reference | | | | 294 | 18776-L27898 | |
| 5q | Reference | | | | 130 | 00797-L00463 | |
| 6p | Reference | | | | 436 | 10731-L11313 | |
| 8q | Reference | | | | 222 | 06746-L27899 | |
| 9p | Reference | | | | 190 | 08067-L19457 | |
| 11p | Reference | | | | 385 | 18677-L30318 | |
| 15q | Reference | | | | 160 | 09787-L10202 | |
| 18q | Reference | | | | 260 | 16433-L27655 | |

SALSA MLPA Probemix P052-D2 Parkinson mix 2

| Chr. position | Target | Exon | Mutation | Distance to next probe | Length (nt) | Probe number | Warning |
|---------------|-----------|---------------------|----------|------------------------|-------------|---------------------|---------|
| 1p36.13 | ATP13A2 | Exon 27 | | 9.4 kb | 448 | 20295-L27934 | |
| 1p36.13 | ATP13A2 | Exon 14 | | | 230 | 11717-L27610 | |
| 4p13 | UCHL1 | Exon 1 | | 0.2 kb | 166 | 03679-L27600 | |
| 4p13 | UCHL1 | Exon 2 | | 0.6 kb | 422 | 21888-L30748 | |
| 4p13 | UCHL1 | Exon 3 | | 3.1 kb | 294 | 20290-L27667 | |
| 4p13 | UCHL1 | Exon 4 | | 1.0 kb | 254 | 20288-SP0953-L28061 | Ж |
| 4p13 | UCHL1 | Exon 5 | | 0.1 kb | 177 | 03681-L03096 | |
| 4p13 | UCHL1 | Exon 6 | | 1.3 kb | 142 | 20285-L27662 | |
| 4p13 | UCHL1 | Exon 7 | | 0.9 kb | 443 | 20294-L27671 | |
| 4p13 | UCHL1 | Exon 8 | | 4.0 kb | 238 | 20287-L27664 | |
| 4p13 | UCHL1 | Exon 9 | | | 203 | 02937-L27602 | |
| 6q26 | PARK2 | Exon 12 | | 10.5 kb | 217 | 06135-L27603 | |
| 6q26 | PARK2 | Exon 11 | | 26.7 kb | 395 | 04614-L27622 | |
| 6q26 | PARK2 | Exon 10 | | 162.2 kb | 350 | 03369-L27619 | |
| 6q26 | PARK2 | Exon 9 | | 20.4 kb | 286 | 20289-L27933 | |
| 6q26 | PARK2 | Exon 8 | | 216.4 kb | 196 | 20286-L27663 | |
| 6q26 | PARK2 | Exon 7 | | 187.5 kb | 161 | 03366-L27599 | |
| 6q26 | PARK2 | Exon 6 | | 80.8 kb | 244 | 03365-L27611 | |
| 6q26 | PARK2 | Exon 5 | | 147.1 kb | 148 | 20257-L27598 | |
| 6q26 | PARK2 | Exon 4 | | 61.5 kb | 343 | 19810-L27618 | |
| 6q26 | PARK2 | Exon 3 | | 180.6 kb | 274 | 05654-L28095 | |
| 6q26 | PARK2 | Exon 2 | | 283.9 kb | 303 | 20291-SP0954-L27668 | Ж |
| 6q26 | PARK2 | Intron 1 | | | 136 | 03204-L02565 | ∅ |
| 7q31.2 | CAV2 | Exon 3 | | 53.0 kb | 477 | 04091-L27626 | |
| 7q31.2 | CAV1 | Exon 3 | | | 404 | 21889-L30747 | |
| 12q12 | LRRK2 | Exon 1 | | 0.4 kb | 379 | 04278-L27621 | |
| 12q12 | LRRK2 | Exon 2 | | 25.9 kb | 486 | 20296-L27936 | |
| 12q12 | LRRK2 | Exon 10 | | 23.2 kb | 429 | 04279-L27624 | |
| 12q12 | LRRK2 | Exon 15 | | 29.3 kb | 466 | 04280-L28024 | |
| 12q12 | LRRK2 | Exon 27 | | 36.3 kb | 281 | 04281-L27614 | |
| 12q12 | LRRK2 | Exon 41 | | 0.1 kb | 190 | 20256-L23585 | |
| 12q12 | LRRK2 | Exon 41 | G2019S | 24.8 kb | 172 | 04574-L27601 | § |
| 12q12 | LRRK2 | Intron 49 (Exon 49) | | | 334 | 04283-L27617 | ∅ |
| 14q22.2 | GCH1 | Exon 6 | | 1.8 kb | 388 | 20292-L27669 | |
| 14q22.2 | GCH1 | Exon 5 | | 1.3 kb | 328 | 03685-L27616 | |
| 14q22.2 | GCH1 | Exon 4 | | 12.6 kb | 369 | 15131-L27620 | + |
| 14q22.2 | GCH1 | Exon 3 | | 5.6 kb | 209 | 04405-L27930 | |
| 14q22.3 | GCH1 | Exon 2 | | 37.7 kb | 319 | 03683-L27615 | |
| 14q22.3 | GCH1 | Upstream (Exon 1) | | | 261 | 04618-L28062 | ∅ |
| 2p | Reference | | | | 500 | 19555-L27674 | |
| 2q | Reference | | | | 154 | 14199-L27215 | |
| 5q | Reference | | | | 130 | 00797-L19287 | |
| 6p | Reference | | | | 359 | 10727-L26803 | |
| 8q | Reference | | | | 224 | 06746-L28025 | |
| 9q | Reference | | | | 415 | 12747-L27779 | |
| 10q | Reference | | | | 310 | 18380-L25673 | |
| 11q | Reference | | | | 183 | 09496-L28060 | |
| 16q | Reference | | | | 268 | 16225-L18478 | |
| 20q | Reference | | | | 460 | 16287-L25505 | |
| 21q | Reference | | | | 494 | 19137-L26747 | |

Probe lengths may vary slightly depending on capillary electrophoresis instrument settings. Please see the most up to date Coffalyser sheet for exact probe lengths obtained at MRC Holland.

The *PINK1*, *ATP13A2*, *PARK7*, *SNCA*, *PARK2*, *LRRK2*, *UCHL1*, *CAV2*, *CAV1*, and *GCH1* exon numbers are derived from the MANE Select transcripts. For more information, see the probe sequences document available on the product page at www.mrcholland.com. Annotations of several probes with targets at the edge of or slightly outside the coding region were altered. The exon numbering from the previous version of this product description is disclosed between brackets.

Chromosomal bands are based on: hg18.

7. Precautions and Warnings

Probe warnings

| | |
|---|--|
| § | These probes will only generate a signal when the mutation is present. |
| - | These probes are flanking probes, included to help determine the extent of a deletion/duplication. |
| Δ | This probe may be sensitive to certain experimental variations. Aberrant results should be treated with caution. |

| | |
|---|--|
| Ж | These probes consist of three parts and have two ligation sites. A low signal of these probes can be due to depurination of the sample DNA, e.g. due to insufficient buffering capacity or a prolonged denaturation time. When this occurs in reference samples, it can look like an increased signal in the test samples. |
| ∅ | These probes target sequences outside of the known coding region. Copy number alterations of only one of these probes are of unknown clinical significance. |

- + The ligation site of this probe is >20 nt away from the nearest exon. For more information, download the probe sequences document available on the product page at www.mrcholland.com.
- ∫ Copy number alterations of this probe are suspected to have clinical significance, but the association is not yet fully established.
- † Copy number alterations of only this probe are of unknown clinical significance.

Probemix-specific precautions

- These products are not known to contain any harmful agents. Based on the concentrations present, none of the ingredients are hazardous as defined by the Hazard Communication Standard. **A Safety Data Sheet (SDS) is not required for these products:** none of the ingredients contain dangerous substances at concentrations requiring distribution of an SDS (as per Regulation (EC) No 1272/2008 [EU-GHS/CLP] and 1907/2006 [REACH] and amendments).
- Sample or technical artefacts may appear as a (mosaic) copy number change of the whole/partial gene. Whole/partial gene deletions or duplications should therefore be confirmed by analysis of an independent DNA sample, to exclude false positive results.
- Small changes (e.g. SNVs, small indels) in the sequence targeted by a probe can cause false positive results, even when >20 nt from the probe ligation site. Sequence changes can reduce the probe signal by preventing ligation of the probe oligonucleotides or by destabilising the binding of a probe oligonucleotide to the sample DNA. Deviations detected by these products should be confirmed, and single-probe deviations always require confirmation. Sequencing of the target region is recommended. Please contact MRC Holland for more information: info@mrcholland.com.
- Copy number alterations of reference probes are unlikely to be related to the condition tested.
- Before testing patient samples, testing of samples from healthy individuals is required to identify suitable reference samples for proper data analysis.

Technique-specific precautions

See the [MLPA General Protocol](#).

8. Limitations

Probemix-specific limitations

- The clinical significance of the following findings is not yet clear/clearly established: multi-exon deletions in *TNFRSF9* may be associated with PD [6].
- The mutation-specific probes can only detect the presence of the mutation and should not be used to determine zygosity.
- Target probes for *PARK7*, *ATP12A2*, *PINK*, *UCHL1*, *GCH1* and *LRRK2* CNVs and the *SNCA* A30P mutation are included to be used for research purposes only and not for diagnostic use.

Technique-specific limitations

See the [MLPA General Protocol](#).

9. References Cited in this IFU

- Schulte, C. and T. Gasser, Genetic basis of Parkinson's disease: inheritance, penetrance, and expression. *Appl Clin Genet*, 2011. 4: p. 67-80.
- Ambroziak, W., et al., *Genomic instability in the PARK2 locus is associated with Parkinson's disease*. *J Appl Genet*, 2015. 56(4): p. 451-461.
- Lesage, S., et al., LRRK2 G2019S as a cause of Parkinson's disease in North African Arabs. *N Engl J Med*, 2006. 354(4): p. 422-3.
- Ozelius, L.J., et al., LRRK2 G2019S as a cause of Parkinson's disease in Ashkenazi Jews. *N Engl J Med*, 2006. 354(4): p. 424-5.
- Ross, O.A., et al., Association of LRRK2 exonic variants with susceptibility to Parkinson's disease: a case-control study. *Lancet Neurol*, 2011. 10(10): p. 898-908.
- Güler, S., et al., Early-Onset Parkinson's Disease: A Novel Deletion Comprising the DJ-1 and TNFRSF9 Genes. *Mov Disord*, 2021. 36(12): p. 2973-2976.

Implemented Changes in the Product Description

Version D2/D2-05 – 24 January 2025 (03S)

- Product description updated to new template.
- Intended purpose was updated; for P051: CNVs for *PARK7*, *ATP13A2*, *PINK* genes and the *SNCA* A30P mutation removed, and for P052: CNVs for *UCHL1*, *LRRK2*, *GCH1*, *ATP13A2*, *CAV1/2* genes removed.
- Salt warning has been removed from ATP13A2 probe 11716-L28013.
- SNV rs566749983 can affect the probe signal. However, the warning for this SNV present in previous product description versions has been replaced by a general warning for small sequence changes, included in section Precautions and Warnings.
- Warning for a ligation site >20nt away from the nearest exon added to probe 15131-L27620.
- Description of probe targets at the edge of or slightly outside the coding region has been adjusted. No change in actual target sites.
- Warning for the target being outside the transcript region added for probes 03690-L27587, 04283-L27617, 04618-L28062, and 03204-L02565.
- Separate warnings for clinical significance added for probes 20271-L27994 and 20224-L27548.
- Probe 03204-L02565 in content table renamed from PACRG exon 1 to PARK2 intron 1 probe.
- Probemixes are now IVDR certified.

MRC Holland, SALSA, MLPA, digitalMLPA, Coffalyser.Net, Coffalyser digitalMLPA, and their logos are trademarks or registered trademarks of MRC Holland BV. All other brands and names herein are the property of their respective owner.